Operating Manual

Laboratory autoclave Steam sterilizer CertoClav Connect

Georg Grinninger-Str. 37 A-4050 Traun / Austria

Content

Introduction4	4
Read Carefully Before Using	4
Proper Use	4
Unwrapping	6
Product contents	6
Place of use ϵ	6
Installation of CertoClav Connect ϵ	6
Assembly and set up CertoTouch controller	7
Turn on CertoTouch	7
Turn on CertoClav Connect pot	7
Choosing the language	7
Checking for updates	7
Connecting CertoTouch with the pot	8
Registering new users	8
Logging in without using the CertoCloud	9
Logging in using the CertoCloud	9
Test run safety instructions	D
Description of the working steps	1
Turning on the CertoTouch1	1
Turning on the CertoClav Connect pot1	1
Connecting the CertoTouch with the the pot	1
Opening lid	2
Pouring in water	2
Loading substance	3
Closing lid	3
Select and run a program	4
Unloading.	5
Interrupting a cycle	6
If power supply is cut off during operation.	6
Description of a program	7
Process status screen	7
Operation of Cycles	7
Sterilization Cycle	7
End of the program	ጸ
Description of the options	R
Flexible media temperature sensor	8
Maintaining of the autoclave	9
Condensate container	9
Cleaning the device housing	9
Real-time-clock battery	9
Gasket	9
Troubleshooting	Ď
Technical data 21	1
Certificates 22	2
Technical data CertoTouch	2
Indication labels	3
Warranty	3

CE Declaration of Conformity

This is to certify, that the **CERTOCLAV** Product Range CertoClav Connect complies with the following regulations of the European Community:

2014/35/EU Low Voltage Directive

Associated Standards:

EN 61010-2-040:2005 Part 2: Particular requirements for sterilizers and washer-disinfectors used to treat medical materials

EN 61010-1:2001 Safety requirements for electrical equipment for measurements, control and laboratory use – Part 1: General requirements

The Certoclav pressure equipment is constructed according Low Voltage Directive.

2014/30/EU Electromagnetic compatibility

Associated Standards:

EN 61326-1:2013 Electrical Equipment for measurements, control and laboratory use - EMC Requirements-Part 1: General requirements

EN 61000-3-2:2006 + A1:2009 + A2:2009 Electromagnetic Compatibility (EMC) — Part 3- 2: Limits Limits for harmonic current emissions (equipment input current equal to or less than 16 A per phase) EN 61000-3-3:2013 Electromagnetic compatibility — Part 3- 3: Limits — Limits for voltage fluctuation and flicker on the public low-voltage system (equipment having an input current equal to or less 16 A per phase) and not subject to conditional connection.

2011/65/EU Restriction of the use of certain hazardous substances (RoHS) Directive

Associated Standards:

EN 50581:2012 Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

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Quality of Production:

CertoClav produces according modern production methods and with high quality materials. Experience for many years and a competent team ensure high quality of products and professional service. All units are 100 % tested on function and quality. Certoclav manufactures in Austria only.

Manufacturer: CERTOCLAV Sterilizer GmbH Georg Grinninger-Str. 37 A-4050 Traun / Austria

(Notes

Traun, 21/11/2016

Mag. Moritz Gruber CEO

Introduction

This document has been carefully compiled. If you have any questions or required information is not covered in the manual, please contact:

CertoClav Sterilizer GmbH Georg Grinninger Str. 37 A-4050 Traun Austria Tel.: +43 (0) 7229 689 Email: support@certoclav.com

Read Carefully Before Using

Please read this manual carefully before using the device to ensure safe and proper use. Please consider that the following instructions insure your own safety and proper use of the device.

Proper Use

The CertoClav Connect small-steam-sterilizer is used to sterilize materials which can withstand temperatures up to 140°C and 2.7 bar of relative pressure. This autoclave is mainly used to sterilize microbiological materials in the laboratory:

- Glass
- Instruments
- Liquid substances and serums
- Media
- Rubber
- Filters
- Synthetic materials

Because of the freely programmable nature od this autoclave it can also be used for material tests and research purposes.

- ▲ Use exclusively the CertoTouch-controller to operate your autoclave.
- ▲ The lid, chamber, gasket and panel are extremely hot immediately after the completion of operation. Do not touch the equipment or you may get burned.
- $\underline{\wedge}$ Keep this manual in a visible, easily accessible location.
- ▲ The autoclave should not be installed in explosive areas
- ▲ Never connect the power cord to a power supply other than one of the rated voltage. Connection to such a power supply can cause fire or electric shock.
- ▲ Only for laboratory use by trained personnel. Do not use this autoclave for medical purposes.
- ▲ This autoclave can be serviced only by our authorised dealers. Only original parts may be installed.
- ▲ Remove the power plug before removing the electronic housing or the bottom.
- ▲ Do not use the autoclave without water. The amount of water needed for sterilization: minimally 0.7 Litre, maximally 2.5 Litre of demineralized soft water.
- ▲ Please do not sterilize tightly closed containers, especially sealed glass bottles without pressure equalisation.
- ▲ Please, transport only with emptied chamber.
- ▲ Do not open the lid under pressure with force.
- ▲ Under no circumstance let the steam in the direction of another person. Risk of a burn. Please use either the provided waste steam container or connect an external waste steam chamber under own responsibility with a hose (inside diameter 10mm) to let out the hot steam.
- ▲ If filled with glass bottles, open the lid only after the temperature of the liquid has dropped under 80°C. Danger of over boil of the liquids or of an explosion of the closed bottles.
- ▲ Under no circumstance put explosive or flammable materials in the autoclave.
- ▲ If the electrical supply is damaged immediately remove the device from the power supply. It can be used again only after a qualified person has repaired the autoclave.
- ▲ Never put objects in the ventilation slot nor pour liquids in the ventilation slot. Danger of electric shocks and short circuits.
- ▲ Manipulating or rebuilding the device or the programs can lead to malfunction.
- ▲ To avoid corrosion in the chamber and the conduit system do no sterilize heavy saline substances.

Unwrapping

- Take the device out of the carton by holding the handles on the lid and lifting it out. Furthermore, remove the CertoTouch controller, power cord and equipment.
- Place the device on stable ground.
- Check if the information on the label coincides with the delivery note.

Product contents

The following items were delivered with the autoclave:

- CertoClav Connect autoclave 18L
- Stainless steel base plate with water level indicator
- Power cord
- CertoTouch control terminal with stand bracket, 4 screws and a power cord
- Operating manual CertoClav Connect
- Test and guarantee certificate

Should one of those parts be missing, please contact your distributor or us directly at: support@certoclav.com or per telephone +43 (0) 7229 689

Should you need additional equipment or a warranty packet, you can order it per telephone, per email or via your distributor or in the online shop at: https://www.certoclav.com/shop/

Place of use

- ▲ Avoid temperatures below 5 ° and beyond 35 °C and air moisture beyond 80 %.
- ▲ The operating terminal is splash-proof. Please guard the autoclave from any splash water.
- ▲ Avoid places where convulsions and vibrations appear and do not use the device outdoors.
- ▲ Take care that the device is not used near steam sensitive laboratory devices (for example microscopes).
- ▲ Avoid places where the ventilation slot might be blocked.

Installation of CertoClav Connect

The voltage values required for this device are to be observed. If your do not know what values the power supply supplies, ask your electrician or at the local electrical utility.

- Turn on the device by connecting the power cord to the power supply and turn the rocker switch on the left side to |.
- If the device was locked, it will now unlock automatically.
- Push and turn the red lever on the right lid handle clockwise.
- Turn the lid anti-clockwise to open the autoclave.
- Take out the baskets and other enclosed equipment.
- Assemble the CertoTouch controller next to the autoclave as described in its operating manual

Assembly and set up CertoTouch controller

- Please mount the stand bracket to CertoTouch-controller with 4 screws. The screws are included in the packaging with the stand bracket.
- If provided, connect the extension cable directly to the CertoTouch controller. There is only one suitable connector.
- Plug the enclosed power supply unit either to the extension cable or directly to the CertoTouch.
- Place the CertoTouch controller only on a levelled, stable surface. A suitable viewing angle can be adjusted by means of the standing bracket.

Turn on CertoTouch

Please press the **U**-button at the front right of CertoTouch for at least 3 seconds.

Turn on CertoClav Connect pot

Turn on by turning the rocker switch to |. If the switch lights green, then the device is turned on. This is necessary to be able to connect the CertoTouch with the CertoClav Connect pot via Bluetooth.

Choosing the language

- 1) Go to the settings menu by clicking the ²-button at the upper right corner of the screen.
- 2) In the Menu settings open the language settings by clicking on the "Language" button.
- 3) In the detail view at the upper right corner of the screen click the current language settings.
- 4) A scrolling list of supported languages can be seen now. Please choose here by clicking the desired one.

Checking for updates

- 1) Please make sure CertoTouch is connected to the Internet.
- 2) Go to the Settings menu by clicking the -button at the top right corner of the screen.
- 3) Please enter the "Device" section by selecting the entry "Device" in the left list.
- 4) In the detailed view on the right, please scroll to "Check for updates" and select this entry.
- 5) A progress bar with the download progress will pop up. Once all updates have been downloaded, they will be installed. If a new software version is available, an installation wizard will guide you through the installation process.

Connecting CertoTouch with the pot

- 1) Please navigate to the Connection Manager. For this purpose, please log out and press the
 - 5

-button at the login window.

- 2) Please, press the ¹-button in order to determine all reachable CertoClav controllers. This may take a few seconds.
- 3) In this list, you can see the serial numbers of all reachable CertoClav controllers.
- 4) Please choose the controller you want to connect with by pressing that serial number in the list. The serial number of each controller can be found on the left or right side of the controller box.
- 5) CertoTouch will connect itself to the selected CertoClav autoclave. The connection is successfully established, if a "connected" status message appears in the selected list sentry. Only after a connection has been successfully established, it is possible to return to the login-window. Please log in there with your email and password in order to use the CertoTouch controller.

Registering new users

- 1) Navigate to the login window. For this, a logged in user must log out first.
- 2) Click on the ¹ button at the top right of the screen.
- 3) Please confirm the popup notification in order to create a new user.

Ð	REGISTER NEW USER	
	Email	
	Password	
	Password	
	REGISTER USER	

- 4) A new window for the user registration will appear.
- 5) Enter the user's email address. Note: The "@" symbol can be reached in the keyboard by pressing the "SYM" (symbols) button.
- 6) Please choose a password with at least 4 characters and retype it.
- 7) Confirm the user registration, by clicking the "Register user" button.

Logging in without using the CertoCloud

1) Make sure that the online mode is deactivated. The online mode can be enabled or disabled in

the settings menu

- 2) Please select the email of the user from the drop-down menu.
- 3) Please type the user password.
- 4) Click the "Login" button.

Logging in using the CertoCloud

- 1) In order to take advantage of the CertoCloud, please create a CertoCloud account first. You can create a CertoCloud account at <u>www.certocloud.com</u>.
- 2) Please connect the CertoTouch with the internet.
- 3) On the CertoTouch, the online mode must be enabled. The online mode can be enabled or

disabled in the settings menu

- Navigate to the login window. For this, a logged in user must log out first.
 Please select the email of the user from the drop-down list. The email must be the same email as the one of the CertoCloud account(<u>www.certocloud.com</u>). If this email address is not available in the dropdown list, refer to the chapter "Register a new user".
- 6) Please enter the CertoCloud password of the user.
- 7) Click the "Login" button.
- 8) Your account is now connected with the CertoCloud. On <u>www.certocloud.com</u> you can access you autoclave, define custom programs and access your protocols.

Test run safety instructions

Please follow the following safety instructions: Please take into consideration all of the instructions and warnings noted on the device.

- ▲ Do not put your face or hands close to the chamber when lifting the lid after a program has finished; steam will gush out of the chamber.
- \triangle Never put your fingers above the steam chamber. The outflowing steam may cause burns.
- △ Open the lid using the handle just a few centimetres first, to let the steam out. Keep safety distance from the device. Afterwards, carefully open the lid fully.
- ▲ Please, keep foreign materials (metals, liquids) from the ventilation, because a current voltage can be induced inside the autoclave, which can cause short circuits or electric shocks. Do not forcibly bend, twist, tie or extend the power cord. Do not place heavy objects on the cord. A damaged cord or exposed wire can cause fire or electric shock.
- ▲ Close the lid after confirming that no foreign matter is adhering to the section contacting the lid gasket. Foreign matter in this section can cause vapor leaks. When using a waste processing bag or other kind of bag and disinfecting, place the bag in the metal mesh holder and then insert it into the chamber. Using the bag "as is" can cause excessive temperatures, pressures, lack-of-water, etc.
- ▲ Please take care not to pinch your hands while closing the lid.
- ▲ The lid, chamber, gasket and panel are extremely hot immediately after completing a program. Do not touch the equipment or you may get burned. Put on heat insulating gloves before removing a substance from the chamber. Do not put hands into the chamber until the steam has been vented.
- ▲ Some time is required for liquids to cool. Be sure to check that the temperature has dropped sufficiently before unloading a liquid from the chamber or burns can result.
- ▲ Do not remove the exhaust bottle before water in the bottle has sufficiently cooled.
- ▲ If any abnormality occurs (e.g. abnormal sounds, smells, smoke), immediately shut off the power. After checking to see that the abnormal condition does not continue, call our authorized distributor in your region.
- ▲ When transporting the device keep the red lever on the handle on "Close" at all times. It prevents the lid from opening on its own.
- ▲ Do not pour anything except demineralized water for sterilizing into the autoclave.

The first test run serves as a functionality control. It should discover the damage that may have occurred during the transport. Furthermore, it demonstrates how to properly handle the device.

Description of the working steps

Turning on the CertoTouch

Please press the **U**-button at the front right of the CertoTouch for at least 3 seconds.

Turning on the CertoClav Connect pot

Turn on the device by turning the switch rocker on |. If the rocker lights green, the device is on. If the temperature in the chamber is under 80°C, the device unlocks itself automatically.

Connecting the CertoTouch with the the pot

1) Please navigate to the connection manager. For this purpose, please log out and press the

-button at the login window.

- 2) Please press the *P*-button in order to determine all reachable CertoClav controllers. This may take a few seconds.
- 3) In this list, you can see the serial numbers of all reachable CertoClav controllers. Please choose the controller you want to connect with, by pressing that serial number in the list. The serial number of each controller can be found at the left or right side of the controller box.
- 4) CertoTouch will connect to the selected CertoClav controller. The connection is successfully established, if a "connected" status message appears on the selected list entry. Only when a connection is successfully established, it's possible to return to the login-window. Please log in there with your email and password in order to use the CertoClav controller.

2016-11-21

Opening lid

Press the red rotary knob on the right handle and turn it clockwise to unlock the lid.

Rotate the lid anti-clockwise to take it off.

Pouring in water

▲ Do not pour anything except demineralized water. Please use soft water since the properties of the water affect the life time of the device.

Please follow the following instructions:

Pour in demineralized soft water in the sterilization chamber until you can see the water through the openings of the base plate. On the base plate is a mark for the minimum amount of water. Needed amount of water: minimally 0.7 Litre, maximally 2.5 Litre.

Loading substance

Please follow the following instructions:

- ▲ Close the lid after confirming that no foreign matter is adhering to the section contacting the lid gasket. Foreign matter in this section can cause vapor leaks. When using a waste processing bag or other kind of bag and disinfecting, place the bag in the metal mesh holder and then insert it into the chamber. Using the bag "as is" can cause excessive temperatures, pressures, lack-of-water, etc.
- ${\ensuremath{\Delta}}$ Be careful not to pinch hands when closing the lid. Always use the base plate.
- ▲ Be sure that the temperature has dropped to at least 50°C in the chamber before the beginning of the next cycle.
- ▲ If a waste disposal bag is used in sterilization, open the bag enough that the bag is not in contact with the inside surface of the chamber. Insufficient sterilization may be caused if the bag is closed during sterilization. When the bag is opened excessively, steam is prevented from circulating in the chamber. This may also result in insufficient sterilization.
- ▲ Do not pile specimen on top of one another. When the chamber is overly packed, steam fails to penetrate all points, resulting in an incomplete sterilization.
- ▲ In sterilizing liquids such as chemicals and media pay attention to the quantity of the liquid in relation to its container. For an Erlenmeyer flask, the amount of chemical should be approx. 3/4 of the capacity of the container; for a test tube, the appropriate quantity of chemical is approx. half of the capacity of the container. Too much chemical may result in overflow from the container during the temperature rising or cooling process.
- ▲ Use container caps that are loose fitting and allow the passage of air. Containers may break if venting is not possible.
- ▲ In the case of dissolution of agar medium, its quantity should be 2 liters or less per container. Two liters or more of agar medium may not be completely dissolved.

Closing lid

Please check that the lid seal is correctly inserted. The signature "Place This Side Facing Down" must be readable. See figure below. Then close the lid in exactly the opposite way, as it was previously opened.

Select and run a program

- 1) Log into your account first. After a user has logged in, the main menu appears.
- 2) Visit the STERILISATION-window by navigation to: MAIN-MENU->STERILIZATION.
- 3) Find the correct program by scrolling the list of available programs.
- 4) By tapping on the start button of the specific program, the user can start the program.
- 5) If activated, please follow the STEP-BY-STEP guide until the program starts. Note: The STEP-BY-STEP guide can be disabled or enabled in the Settings->Sterilisation->STEP-BY-STEP-GUIDE section.

For further instructions on how to use the CertoTouch please refer to the provided operating manual.

Unloading

In order to open the device and take out the sterilized goods must the temperature in the chamber drop to under 80°C. This serves as a safety measure to prevent burns. As soon as the cycle ends, a message "Please open the lid" appears on the screen. The program counts as successfully finished only after opening of the lid.

- ▲ Keep the face and hands away from the chamber when lifting the lid after operations are complete; steam will gush out of the opening of the chamber.
- ▲ After operation is finished, the lid, chamber, packing and panel are hot. To protect yourself from burns, do not touch with bare hands. Put on heat insulation gloves before taking a substance out of the chamber. Do not put your hands into the chamber until the steam has been vented.
- ▲ It takes a lot of time for the liquid to cool. Be sure to check that the temperature has dropped sufficiently before unloading the liquid from the chamber, or you may get burned.

To open the lid please follow the instructions:

Press the red rotary knob on the right handle and turn it clockwise to unlock the lid.

Rotate the lid anti-clockwise to take it off.

Interrupting a cycle

- ▲ If the cycle is interrupted during the sterilization of liquids, the lid is electrically unlocked after the temperature in the chamber has dropped to under 80°C.
- ▲ In case of power loss, the lid locks itself.

If power supply is cut off during operation

If the power supply is cut off due to power failure or the like, operation is interrupted. When power supply is restored, the equipment will be back in the standby state (state before operation). In this case, repeat operations from the beginning. The autoclave can not be opened when turned off.

Description of a program

Process status screen

If a program is running, the MONITOR window is visible. In the monitor window, all important information about the current cycle is displayed. The following information can be found (additionally marked in the screenshot below).

Name of the running program.

Status of the running program. This can be "Running", "Finished", "Prepare to start", "Not running".

Description of the running program.

Current temperature in [°C] and pressure in [kPa], measured in the chamber.

Graphical overview.

Operation of Cycles

3 4

If activated in the program, the autoclave will first be aired.

Any air remaining in the chamber makes the temperature distribution in the chamber uneven. This hinders temperature rise (sterilization). To cope with this, a microcomputer-controlled automatic exhaust valve vents out virtually 100 % of the air in 10 minutes at 115°C. In the chamber remains water steam. At the end the valve closes.

Sterilization Cycle

A constant temperature (pressure) is maintained during the set sterilization time period. If the temperature in the chamber drops 1°C or more or rises 2°C or more from the set value due to any trouble, the program is interrupted.

▲ In sterilization of petri dishes or empty containers, the air remaining in the container expands and may increase the pressure extraordinarily within the chamber. If the pressure in the chamber exceeds the saturated steam pressure, the automatic exhaust valve opens and discharges the air in the chamber into the exhaust bottle.

End of the program

After all the cycles have been finished, a message "Please, open the lid" comes forth. When the lid is opened, the program is finished.

▲ If the touch screen is not touched for 30 minutes, the display becomes black. By touching the touch screen will the display turn on again.

Description of the options

Flexible media temperature sensor

In general

The flexible media temperature sensor measures the temperature in the media that need to be sterilized. Because of the delay of boiling, for example by sterilization of liquids, it is reasonable to use a media temperature sensor. Additionally, you can control the exact temperatures, which the media was exposed to, after the sterilization cycle.

- ▲ Do not forcibly bend, twist, tie or extend the cable of the temperature sensor. Do not place heavy objects on the cord. A damaged cord or exposed wire can cause malfunction.
- ▲ When sterilizing plastic petri dishes, please, take care that the media temperature sensor does not contact any part of the petri dishes. Because the plastic sclerotizes when the temperature drops, it becomes impossible to separate those parts from the temperature sensor.

Use

After loading the chamber, set the media sensor in the middle of the media you want to sterilize. The media temperature during the program always has precedence over the chamber temperature. For example, the sterilization duration begins only after the media has reached the sterilization temperature.

Deactivating of the media sensor

In the program settings (see CertoTouch operating manual) it can be specified whether the media sensor should be ignored or used.

Maintaining of the autoclave

Condensate container

During every sterilizations process the water gathers in the condensate container. Empty the condensate container after each cycle. Check before each new cycle, if the condensate container has been emptied.

- ▲ Do not empty the condensate container before it has cooled sufficiently.
- ▲ The steam exhaust valve closes itself if the condensate container was not installed.

Cleaning the device housing

▲ Please, do not use any aggressive detergents to clean the device housing.

Real-time-clock battery

The devices intern battery supplies a real-time-clock with power. As soon as the battery dies, there will be a warning showed on the screen that must be confirmed. The average life expectancy of the battery is 5 years. The autoclave can also be operated with a dead battery. Without the real-time-clock it may come to inaccuracies up to a couple of minutes at longer user-defined programs. To change the battery please contact us. When your autoclave is not connected to the network or the internet there is a possibility of fluctuations in time and date.

Gasket

A brittle gasket could lead to leakage of steam. The average life expectancy of a gasket is 4 years. When the gasket needs to be replaced or you wish to buy a new one, please contact us.

Troubleshooting

Should any malfunction occur in the autoclave, the error detection circuit will be triggered to assure safety. Once the circuit is activated, an error number appears on the digital display. The alarm will stop only after the error has been confirmed.

Error code	Error description	Next step	
Error 1	Faulty ventilator	Replace the ventilator	
Error 2	The lid is not closed	Please contact the service-team.	
		Please, close the lid prior to staring a program. Follow	
		the step-by-step video.	
Error 3	Too much water in	Please, let some water out of the chamber. There	
	the chamber	should be maximally 2,5 Litre water in the chamber	
		during the cycle. If the error does not disappear, please	
		contact the service-team.	
Error 4	Too little water in	Please, fill the chamber with minimally 0,7 Litre water in	
	the chamber	the chamber before the start of a program. If the error	
		does not disappear, please contact the service-team.	
Error 5	The lid cannot be	Please, turn off the device and try once more. If the error	
	locked	does not disappear, please contact the service-team.	
Error 6	Sterilization heater	Temperature in the chamber cannot be achieved. Either	
	trouble alarm	the autoclave is empty of water or the heater is faulty.	
		Please, contact the service-team.	
Error 7	Over – pressure	Please, contact the service-team.	
	alarm		
Error 8	Over – temperature	Please, contact the service-team.	
	alarm	Possible defects: temperature sensor is faulty or the	
		heater malfunctions.	
Error 9	Pressure out of	Unexpected pressure fluctuations in the chamber.	
	permissible range	Please, do not put any tightly closed containers in the	
		chamber.	
Error 10	Temperature in the	Please, contact the service-team.	
	chamber out of	Temperature fluctuations during the cycle.	
	permissible range	Possible causes: heater or temperature sensor broken.	
Error 11	Temperature of the	Please, contact the service-team.	
	substance out of	Temperature fluctuations during the cycle.	
	permissible range	Possible causes: heater defect, temperature sensor	
		defect.	
Error 13	Temperature band	Please, contact the service-team.	
	not adhered to	Possible combined error of Error 10 and Error 11	
Error 14	Cycle aborted	If a fatal error during the sterilization occurs, and the	
	because of a fatal	program must be aborted, then will an additional error	
	error.	be showed.	
Error 15	Cycle aborted	No action needed. The program was manually	
	because of user's	cancelled. Please, wait until the temperature in the	
—	interference.	chamber has dropped to under 80°C.	
Error 16	Tablet-PC and	Please set the Touch controller closer to the autoclave	
	autoclave	and renew the connection. You should not move the	
	connection	I ouch controller further that 6m from the autoclave.	
	terminated		

Technical data

Measures and weight			
Material	Stainless steal	Certificates	CE
	304L		
Dimensions			
Chamber Volume	18 Litre	Device height	60 cm
Chamber diameter	26 cm	Device width	44 cm
Chamber height	35 cm	Device depth	40 cm
Chamber usable height	30 cm	Device weight	15 kg
Operating data			
Max. operating	140 °C	Power supply	230V, 50/60 Hz
temperature			
Max. operating pressure	2,7 Bar (39 psi)	Power	1300 W
Test pressure	4,1 Bar (59 psi)		

Heating duration

Amount of load	Heating from 25 °C to 121 °C
Little	15 – 20 Minutes
A lot	25 – 40 Minutes

Cooling duration

Amount of load	Cooling duration at faster steam exhaust from 121 °C to 80°C	Cooling duration at slower steam exhaust from 121 °C to 80°C
Little	20 – 35 minutes	40 – 60 minutes
A lot	40 – 60 minutes	80 – 120 minutes

Correlations table

Tin℃	p in bar (kgf/cm	T in ℃	p in bar (kgf/cm²
100	0,03	118	0,90
101	0,07	119	0,96
102	0,11	120	1,02
103	0,15	121	1,09
104	0,19	122	1,15
105	0,23	123	1,22
106	0,28	124	1,29
107	0,32	125	1,37
108	0,37	126	1,44
109	0,41	127	1,52
110	0,46	128	1,59
111	0,51	129	1,67
112	0,56	130	1,75
113	0,61	131	1,84
114	0,67	132	1,92
115	0,72	133	2,01
116	0,78	134	2,10
117	0,84	135	2,19

Certificates

CE certified, Please see chapter EC declaration of conformity.

Technical data CertoTouch

Display screen	7" LED background light	
Resolution	480 pixels	
Touch panel	capacitive	
Brightness	450cd/m ²	
Aspect ratio	16:9	
CPU	Freescale iMX.53 x 800MHz / 1.0GHz	
Core	ARM Cortex-A8	
Memory card connectors	Micro-SD card und USB-Stick optional	
LAN (cable)	Ethernet 10/100M (RJ45) optional	
W-LAN (without cable)	IEEE 802.11	
Bluetooth	Bluetooth 2.1 Host	
Other Interfaces	CertoClav barcode scanner port	
CertoClav label printer port		
Audio Earphone slot		
	Speaker	
Speaker	1W	
Power consumption	9W	
Working temperature	-10°C to 65°C	
Storage temperature	-25°C to 80°C	
Power supply	DC 9V to 36V	
Dimensions (LWD)	220x132x36.5mm	
Weight	740g	
Degree of protection	All-round splash water protection	

	Directive 2002/96/EC (WEEE): the product should not be disposed of with your other household waste. It is your responsibility to dispose of your waste equipment separately from the municipal waste stream.		
\bigwedge	WARNING: Read the operating manual before usage		
	WARNING: Hot area		

Indication labels

Warranty

This autoclave was quality, safety and functionality checked by our company. We provide a 2-year warranty on all parts, that are not consumables, from the day of the delivery.

Please keep the separate warranty card for future reference.

An optional 5-year extension to the warranty can be purchased.

No warranty in case of inappropriate use.

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CE certified according to: Low voltage directive 2006/95/EC EN 61010-1:2001. EN 61010-2-40:2005 (IEC 1010-1:2001. IEC 1010-2-040:2005) ZEK 01.4-08/11.11

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Display screen	7" LED background light
Resolution	480 pixels
Touch panel	capacitive
Brightness	450cd/m ²
Aspect ratio	16:9
CPU	Freescale iMX.53 x 800MHz / 1.0GHz
Core	ARM Cortex-A8
Memory card connectors	Micro-SD card und USB-Stick optional
LAN (cable)	Ethernet 10/100M (RJ45) optional
W-LAN (without cable)	IEEE 802.11
Bluetooth	Bluetooth 2.1 Host
Other Interfaces	CertoClav barcode scanner port
	CertoClav label printer port
Audio	Earphone slot
	Speaker
Speaker	1W
Power consumption	9W
Working temperature	-10°C to 65°C
Storage temperature	-25°C to 80°C
Power supply	DC 9V to 36V
Dimensions (LWD)	220x132x36.5mm
Weight	740g
Degree of protection	All-round splash water protection

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