



## SPECIFICATIONS


|                       | FPV 01 (Low Viscosity)  | FPV 02 (High Viscosity)  |
|-----------------------|---|--|
| Measurement Range     | No. 4 rotor: 2 to 33 mPa·s<br>No. 5 rotor: 15 to 150 mPa·s<br>No. 3 rotor: 50 to 300 mPa·s  | No. 3 rotor: 0.3 to 13 dPa·s (with No. 3 cup)<br>No. 1 rotor: 3 to 150 dPa·s (with JIS 300 mL beaker)<br>No. 2 rotor: 100 to 4000 dPa·s (with JIS 300 mL beaker)   |
| Sample Fluid Capacity | Approx. 460 mL (with Cup A or Cup B)  | No. 1 and No. 2 rotor (with JIS 300 mL beaker) approx. 300 mL<br>No. 3 rotor (with No. 3 cup) approx. 170 mL<br>Clearance between rotor end and cup bottom: about 15 mm.   |
| Measurement Accuracy  | Within ±5% of maximum measurement range for each rotor  | ±10% of indicated value, reproducibility ±5%   |
| Rotor Speed           | 62.5 rpm  |  |
| Power Supply          | IEC LR6 (size AA) alkaline batteries, nickel-hydrate rechargeable batteries, AC adapter 0.5A  |  |
| Dimensions and Weight | 175 (H) x 77 (W) x 40 (D) mm (without protruding parts). Approx. 260 g (without batteries)  |  |
| Supplied Accessories  | No. 3 rotor (dia 45x 47 x 160 mm) SUS304 1<br>No. 4 rotor (dia. 78 x 46 x 159 mm) A1050 (alumite) 1<br>No. 5 rotor (dia. 61.2 x 36 x 149 mm) A1050 (alumite) 1<br>Cup A (dia. 92 x 76 mm, without hole) A1050 (alumite) 1<br>Cup B (dia. 92 x 76 mm, with hole) A1050 (alumite) 1<br>IEC LR6 (size AA) alkaline batteries 4 | No. 1 rotor (dia. 24 x 53 x 166 mm) SUS304 1<br>No. 2 rotor (dia. 15 x 1 x 113 mm) SUS304 1<br>No. 3 rotor (dia. 45 x 47 x 160 mm) SUS304 1<br>No. 3 Cup (dia. 45 x 47 x 160 mm) SUS304 1<br>Extension rod (900mm, 300x3) SUS304 1<br>IEC LR6 (size AA) alkaline batteries 4 |

1 JIS R 3503 : 1994 78x103

### Options

| Product Name | Product Number |
|--------------|----------------|
| Stand        |                |
| AC Adapter   |                |

## ROTORS AND CUPS (UNIT: mm)



## VISCOTESTER MEASUREMENT EXAMPLES (FOR REFERENCE)

| Product Type                 | Viscosity  | Viscotester | Rotor |
|------------------------------|------------|-------------|-------|
| <b>Newtonian Fluids:</b>     |            |             |       |
| Milk                         | 2.6 mPa·s  | FPV 01      | No. 4 |
| Soy Sauce                    | 5 mPa·s    | FPV 01      | No. 4 |
| Lactic Fermented Beverage    | 28 mPa·s   | FPV 01      | No. 5 |
| Olive Oil                    | 71 mPa·s   | FPV 01      | No. 5 |
| Castor Oil                   | 6 dPa·s    | FPV 02      | No. 3 |
| Starch Syrup                 | 1000 dPa·s | FPV 02      | No. 2 |
| <b>Non-Newtonian Fluids:</b> |            |             |       |
| Tomato Juice                 | 230 mPa·s  | FPV 01      | No. 3 |
| Condensed Milk               | 16 dPa·s   | FPV 02      | No. 1 |
| Chocolate Syrup              | 25 dPa·s   | FPV 02      | No. 1 |
| Tomato Ketchup               | 43 dPa·s   | FPV 02      | No. 1 |
| Pure Honey                   | 76 dPa·s   | FPV 02      | No. 1 |
| Toothpaste                   | 320 dPa·s  | FPV 02      | No. 2 |
| Starch Paste                 | 310 dPa·s  | FPV 02      | No. 2 |

|                                      | FPV 01         | FPV 02         |
|--------------------------------------|----------------|----------------|
| Cup A                                | approx. 460 mL | —              |
| No. 3 Cup                            | —              | approx. 170 mL |
| Commercially available 300 mL beaker | —              | approx. 350 mL |

Note: For certain fluids, readings may differ slightly from other viscometers, depending on properties of target fluids, mechanical factors, as well as specific gravity, rotor speed, and other aspects.

$$1 \text{ cP} = \frac{1}{1000} \text{ Pa} \cdot \text{s} = 1 \text{ mPa} \cdot \text{s}$$

$$1 \text{ P} = \frac{1}{10} \text{ Pa} \cdot \text{s} = 1 \text{ dPa} \cdot \text{s}$$