MIKRO 120

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Gartenstraße 100 • D-78532 Tuttlingen • Germany

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Contents

1	Intended application.....	13
2	Notes on safety.....	13
3	Information symbol on the centrifuge.....	14
4	Delivery checklist.....	14
5	Space requirement	14
6	Connection to the mains	14
7	Commissioning	14
8	Opening and closing the lid	14
8.1	Opening the lid	14
8.2	Closing the lid	14
9	Installation and removal of the rotor.....	15
10	Arming the rotor	15
11	Control panel.....	15
11.1	Speed area	15
11.2	Rotation area	15
11.3	Time area.....	15
11.4	Key area	15
12	Adjustable parameters	16
13	Preselecting centrifuging parameters or changing them during operation	16
13.1	Speed	16
13.2	- of denser substances.....	16
13.3	Run time	16
13.3.1	Continuous operation.....	17
13.3.2	Pulsed operation	17
13.4	Brake adjustment.....	17
13.5	Start centrifuging run.....	17
13.6	End centrifuging run.....	17
14	Relative centrifugal force (RCF).....	17
15	Emergency release	18
16	Care / maintenance.....	18
17	Faults	19
18	Change mains input fuse	19
19	Acceptance of the centrifuges for repair.....	19
20	Anhang / Appendix.....	36
20.1	Technische Daten / Technical specification	36
20.2	Rotoren und Zubehör / Rotors and accessories.....	37

1 Intended application

The centrifuge is used for separating substances or mixtures with a density of up to max. 1.2 kg/dm³.

2 Notes on safety

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner.
These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate.
 - However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge should be installed on a good, stable base.
- Ensure that no persons and / or dangerous materials are located within a safety zone of 300 mm around the equipment when a centrifuge is running.
- The centrifuge must not be moved or knocked during operation.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either run hot in the cold room for 30 minutes, or warm up for 3 hours in the warm room, before connecting to the mains.
- Load centrifuge rotor evenly. All positions on rotor must be filled.
- Do not fill centrifuge containers inside the centrifuge.
- Centrifuge containers must not be filled beyond the capacity specified by the manufacturer.
- Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2)
- Centrifuge containers may only be centrifuged with accessories (reducing adapters, frames, suspensions, etc.) authorised by the manufacturer (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge must not be operated in areas subject to danger of explosions.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
Without additional proceedings (like an additional bioseal between bucket and lid of bucket or angle rotor with a special bioseal between rotor and lid) a centrifuge is not a biosafety system in accordance to the regulation EN 61010-2-20. In the case of material belonging to risk group II (see the World Health Organisation's "Laboratory Biosafety Manual") they should employ a biosafety system. Under this system small drips and aerosols are prevented from escaping by a bioseal (packing ring) located between the hanger and the lid. Centrifuge containers with special screw caps, as obtainable through trade suppliers, can also be used for hazardous substances.
In the case of materials from the higher risk groups greater safety provision is required than the arrangements described above. In a biosafety system, centrifuge containers with special screw caps must be used.
- For further details of available biosafety systems see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Any rotors, hangers or accessories showing clear signs of corrosion or mechanical defects must not be used for centrifuging.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Hettich company are allowed to be utilised.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- This centrifuge is classified in Germany as a Group 3 device according to the *Medizinische Geräteverordnung MedGV* (the regulations on medical equipment).
- It conforms to safety regulations based on:
 - IEC 1010-1/-2, DIN - EN61010 Parts 1and 2
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of IEC stipulations,
 - prescribed tests to BGV A1, BGR 261 are carried out by an expert.
- With centrifuges for robotic use please pay attention to the notes of the key operated switch.

No claim under guarantee will be considered by the manufacturer unless the above instructions have been adhered to.

3 Information symbol on the centrifuge

Before the initial operation of your centrifuge you should read and pay attention to the operating instructions.

4 Delivery checklist

The following items and accessories are delivered with the centrifuge:

- 1 Connecting cable
- 2 Fuses
- 1 Notes on moving the equipment safely
- 1 Operating instructions
- 1 Hex. pin driver for transport protection

The rotor(s) and associated accessories are included in the delivery in the quantity.

5 Space requirement

- According to the laboratory instrument standards EN 61010-2-20 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system. This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- The necessary space requirement can be found under Dimensions in the Technical data chapter.
- The centrifuge must be set up in a suitable place, so that it is stable. During set up the required safety area of 300 mm around the centrifuge, in accordance with IEC 1010-2-2, must be observed.

Persons and hazardous materials must not be located in the safety area whilst the centrifuge is in operation.

- Do not place any object in front of the ventiduct.
 - Keep a ventilation area of 300 mm around the ventiduct.

6 Connection to the mains

- Check whether the supply voltage, supply frequency and on-site mains fuse agree with the specification on the nameplate. The nameplate is located on the back of the centrifuge.
- Make sure that the mains switch is in the "0" position.
- The centrifuge must be connected to a standard mains socket using the power supply cable provided.

7 Commissioning

- Check that the centrifuge has been properly set up and that the electrical connections are correct (see Connection to the Mains and Space Requirement).
- Switch the mains switch "ON", switch position "I". After a short time the control panel will switch on and the last set parameters will appear in the displays.
- When the  symbol has lit up, the lid can be opened.
- Open the lid and remove the transport protection (see Transport Protection information sheet).

8 Opening and closing the lid

8.1 Opening the lid

- When the  symbol has lit up, the lid can be unlocked and opened.
- Unlock the lid by pressing the  key and swivel the lid upwards.
- The  symbol appears.

The lid can only be opened if the centrifuge is switched on and the rotor is stationary. If this is not possible, see chapter "Emergency release".


8.2 Closing the lid

- Close the lid by swinging the lid down and pressing lightly. If the symbol  lights up in the rotation display the lid has been correctly locked.

Do not bang the lid shut.

9 Installation and removal of the rotor

- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The driver areas (D) of the motor shaft must be located in the groove (B) of the rotor. The alignment of the groove is labelled on the rotor.
- Tighten up the clamping nut of the rotor by turning it clockwise by hand.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tensioning nut by turning anti-clockwise. Turn the tensioning nut until the rotor is able to be lifted from the motor shaft.


10 Arming the rotor

- Always fill the centrifuge containers outside the centrifuge.
- Ensure, by eye, that the containers are filled to a uniform level.
- Load positions located opposite one another equally.
- Ensure that the rotor is correctly positioned and secured.
- The weight of the permissible filling quantity is specified on each rotor. This weight may not be exceeded.

11 Control panel

The control panel is divided into four areas.

- Speed
- Rotation display
- Time
- Keys

11.1 Speed area

The speed area contains the speed display which indicates the preselected speed or the ACTUAL speed during operation.

The speed can be preselected or changed during operation with the ▲ and ▼ arrow keys. At the end of operation, the display flashes in seconds until the lid is opened or a key is activated.

11.2 Rotation area

The rotation area contains the rotation display, which illuminates and rotates anti-clockwise whilst the rotor is turning.

When the ⏸ symbol lights up, the rotor is stationary and the lid can be opened.

When the lid is unlocked, the ⏸ symbol lights up.

11.3 Time area

This area contains the running time display and the ▲ and ▼ arrows, with which the switch-on time can be preselected or changed during operation.

When stationary, the preselected switch-on time is displayed, as well as the remaining run time during operation. The decimal point flashes each second during operation.

At the end of operation, the display flashes each second until the lid is opened or a key is activated.

11.4 Key area

Start key

To start a run with the preselected parameters.

Stop key 	The STOP key allows a centrifuging run to be stopped at any time.
Pulse key 	Brief centrifuging whilst the key is pressed down. The run time is displayed in seconds in the run time display. After 60 seconds, the run time is displayed in minutes.
Lid key 	After standstill the lid can be unlocked by pressing the corresponding key.

12 Adjustable parameters

Speed	<ul style="list-style-type: none"> • min. speed 500 rpm • max. speed 14000 rpm • in steps of 100 rpm Adjustable up to max. speed of the rotor used.
Time " _ " 	<ul style="list-style-type: none"> • min. preset time 1 min • max. preset time 99 min • in steps of 1 min or <ul style="list-style-type: none"> • continuous operation or <ul style="list-style-type: none"> • pulse (short-time)
Braking force	<ul style="list-style-type: none"> • full braking effect • low braking effect

13 Preselecting centrifuging parameters or changing them during operation.

13.1 Speed

Each time the  and  arrow keys are activated, the speed is changed by 100 rpm. If an arrow key is kept pressed down for longer, the speed is changed at an increasing rate.

If the speed is changed during operation, this becomes effective immediately, i.e. the speed is adjusted.

13.2 - of denser substances

The rotors are designed to centrifuge substances up to a maximum mean homogenous density of 1.2 kg/dm³ when rotating at the stated speed.

Denser substances must be centrifuged at lower speed.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density}}} \times \text{Rated speed}$$

e.g.: RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2}{1.6}} \times 4000 = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

13.3 Run time

Each time the  and  arrow keys are activated, the run time is changed by 1 minute. If an arrow key is kept pressed down for longer, the run time is changed in minute steps, at an increasing rate.

If the run time is changed during operation, this becomes effective immediately, i.e. the remaining run time is adjusted.

If an arrow key is activated during operation, the run time or speed is adjusted accordingly.

13.3.1 Continuous operation

The centrifuge can be operated in continuous operation. Continuous operation is set by pressing the arrow key until “–” appears in the run time display. Operation is only stopped by activating the key.

13.3.2 Pulsed operation


For brief centrifuging. The rotor turns at the preselected speed as long as the key is pressed down.

13.4 Brake adjustment

With this centrifuge, the braking effect can be set to normal or low.

This can be set before a run, as follows:

1. Switch off the mains switch
2. Simultaneously press the pulse key and the arrow key in the speed area and switch on the mains switch. Release the keys again.
3. If required press the arrow button in the speed field until the depicted display appears. In the speed field the machine version as set by the factory and in the time field the brake setting.

4. Set the required value, “1” or “0”, in the time area with the and arrow keys.
5. Confirm the set value with the key.
See Rotors and Accessories in the Braking Times chapter.

13.5 Start centrifuging run

After setting the centrifuging parameters, start the centrifuge by pressing the key.

13.6 End centrifuging run

A centrifuging run is ended by pressing the key or after expiry of the set run time. Once the rotor has come to a stop the lid lock is briefly triggered and the display flashes. The lid can now be opened.

14 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$RCF = \left(\frac{RPM}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad RPM = \sqrt{\frac{RCF}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

15 Emergency release


If loss of power or a centrifuge fault occurs while the centrifuge is running, the lid remains locked.

To release the lid in an emergency, unplug the centrifuge from the power.
Wait for the rotor to stop turning before opening the lid.
To release the lid in an emergency use the plastic release pin, attached at the bottom of the centrifuge.
Important: Use only the attached pin to release the lid.

The release pin is fixed with a clamp at the bottom of the centrifuge and tied to it with a plastic cord.

1. Pull out the release pin and uncoil the cord, see drawing.
2. Insert the release pin horizontally into the hole located in the middle of the front panel, see drawing.
3. Push the release pin to the right, then up to a 45° position. Pull the lid upwards at the same time to unlock, see drawing.
4. Pull the release pin out and recoil the cord around the clamp. Refasten the release pin as shown below.

16 Care / maintenance

Pull the mains plug before cleaning.
Before applying any cleaning or disinfecting procedure other than those recommended by the manufacturer, the user should contact the manufacturer to make sure that the planned process will not damage the equipment.

- The centrifuge should be cleaned regularly for reasons of hygiene, and if necessary should also be cleaned with soap or a mild cleaning agent.
- Any adherent impurities should be removed as they can cause corrosion.
- Humidity in the air or centrifuge containers with no hermetic seal can lead to condensation. The centrifuge chamber should therefore be cleaned regularly with a cloth or similar.
- In the case of glass breakage, the fragments of glass along with any spilt centrifuge product should be removed carefully from the centrifuge chamber, the containers or container drill holes.

After a glass breakage the rubber inserts for the containers must be replaced because any residual glass fragments in these inserts can cause further glass breakage.

- If any infectious material should find its way into the centrifuge chamber it should be disinfected immediately.
- The cleaning and disinfectant agents used should be in a pH range between 5 and 8 (Helipur H plus N, B. Braun Melsungen). Alkaline cleaning agents with pH value > 8 should be avoided.
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

17 Faults

If the fault cannot be rectified according to the faults table, Hettich customer services must be informed.

Please state the type of centrifuge and the factory serial number. Both values are visible on the centrifuge type plate.


MAINS RESET: - Mains switch OFF for longer than 10 secs.
- Mains switch ON.

Fault	Display	Cause of fault	Remedy
No display	---	No voltage Triggering of the excess current cut-out.	- Check distribution voltage. - Check fuses. - Mains switch ON.
System reset	- 2 -	Failure of power supply during a run Unbraked runout.	- When stationary, open lid and press START key.
Balance error	- 3 -	Balance error on the motor axis due to weight differences in the rotor loading.	- Open lid when rotor is stationary. - Eliminate balance error.
Kommunication	- 4 -	Fault in control unit or power unit. Unbraked runout.	- Perform mains reset when stationary.
Overload	- 5 -	Motor or motor control defective.	- Perform mains reset when stationary.
Overvoltage Undervoltage	- 6 - - 8 -	Supply voltage outside tolerance (see Technical Data) Unbraked runout.	- Perform mains reset when stationary. - Check supply voltage.
Excess temperature	- 9 -	Excess temperature switch in motor has triggered. Unbraked runout.	- When rotor is stationary, open lid using emergency unlocking (see Emergency Unlocking chapter). - Allow motor to cool down.
Controller watchdog	- C -	Fault in control unit Unbraked runout.	- Perform mains reset when stationary.
Lid error	- d -	Unbraked rundown, after standstill lid release	- Perform mains reset after standstill.
Short circuit	- E -	Short circuit in control unit / power unit.	- Perform mains reset when stationary.
Configuration	- F -	Incorrect machine version	- Inform Customer Service.

18 Change mains input fuse

Switch off the mains switch and separate the centrifuge from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuse.
- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

19 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

You will find further information on our homepage: <http://www.HettichLab.com>

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No part of this publication may be reproduced without the written prior permission of the copyright owner.