**MIKRO 200  
MIKRO 200 R**

(DE) <b>Bedienungsanleitung .....</b>	<b>5</b>
(EN) <b>Operating Instructions .....</b>	<b>17</b>
(FR) <b>Mode d'emploi .....</b>	<b>29</b>
(IT) <b>Istruzioni per l'uso .....</b>	<b>41</b>


Fig. 1


Fig. 2 MIKRO 200


Fig. 3 MIKRO 200 R

**EG-Konformitätserklärung  
EC Conformity Declaration  
Déclaration de conformité CE  
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

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## Contents

1	Use according to specification .....	19
2	Residual risks .....	19
3	Notes on safety.....	19
4	Explanations of symbols and references .....	20
4.1	Symbol used on the machine .....	20
4.2	Symbols used in the operating instructions .....	20
5	Delivery checklist.....	20
6	Unpacking the centrifuge .....	20
7	Initial operation .....	21
8	Opening and closing the lid.....	21
8.1	Opening the lid .....	21
8.2	Closing the lid .....	21
9	Installation and removal of the rotor.....	21
10	Loading the rotor .....	22
11	Aerosol tight sealing of angle rotors.....	22
12	Control and display elements .....	22
12.1	Control knob.....	22
12.2	Control panel pushbuttons (keys) .....	22
12.3	Adjustment possibilities.....	23
13	Programming .....	23
13.1	Programme input/alteration.....	23
13.2	Programme recall.....	23
14	Centrifugation .....	23
14.1	Centrifugation with pre-set time .....	24
14.2	Continuous run.....	24
14.3	Short-term centrifugation.....	24
15	Emergency Stop .....	24
16	Acoustic Signal .....	24
17	Recall hours of operation.....	24
18	Cooling (only in centrifuges with cooling).....	25
18.1	Standby-cooling .....	25
18.2	Pre-cooling the rotor.....	25
19	Relative centrifugal force (RCF) .....	25
20	Centrifugation of materials with higher density .....	26
21	Emergency release.....	26
22	Maintenance and servicing .....	26
22.1	Centrifuge .....	26
22.2	Rotors and Attachments.....	27
22.3	Autoclaving .....	27
22.4	Centrifuge containers .....	27
23	Faults.....	28
24	Change mains input fuse .....	28

---

25	Acceptance of the centrifuges for repair .....	28
26	Anhang / Appendix .....	53
26.1	Technische Daten / Technical specification .....	53
26.2	Rotoren und Zubehör / Rotors and accessories.....	54

## 1 Use according to specification

The machine on hand is a centrifuge that is exclusively conceived for the separation of substances or respectively substance mixtures with a density of max. 1.2 kg/dm<sup>3</sup> and which is therefore only designated for this application. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

## 2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. However during your application of the device hazards towards life and limb of the user or third parties or respectively impairments to the device or other material assets can arise. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

## 3 Notes on safety

- Before the initial operation of your centrifuge you should read and pay attention to the operating instructions.  
**Only personnel that has read and understood the operating instructions are allowed to operate the device.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner.  
These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate.
  - However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- Ensure that no persons and / or dangerous materials are located within a safety zone of 300 mm around the equipment when a centrifuge is running.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either run hot in the cold room for 30 minutes, or warm up for 3 hours in the warm room, before connecting to the mains.
- Load centrifuge rotor evenly. All positions on rotor must be filled.
- Do not fill centrifuge containers inside the centrifuge.
- Centrifuge containers must not be filled beyond the capacity specified by the manufacturer.
- Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2)
- Centrifuge containers may only be centrifuged with accessories (reducing adapters, frames, suspensions, etc.) authorised by the manufacturer (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge must not be operated in areas subject to danger of explosions.
- The centrifuge must not be used with:
  - inflammable or explosive materials
  - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.  
Without additional proceedings (like an additional bioseal between bucket and lid of bucket or angle rotor with a special bioseal between rotor and lid) a centrifuge is not a biosafety system in accordance to the regulation EN 61010-2-020. In the case of material belonging to risk group II (see the World Health Organisation's "Laboratory Biosafety Manual") they should employ a biosafety system. Under this system small drips and aerosols are prevented from escaping by a bioseal (packing ring) located between the hanger and the lid. Centrifuge containers with special screw caps, as obtainable through trade suppliers, can also be used for hazardous substances.  
In the case of materials from the higher risk groups greater safety provision is required than the arrangements described above. In a biosafety system, centrifuge containers with special screw caps must be used.
- For further details of available biosafety systems see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.

- Any rotors, hangers or accessories showing clear signs of corrosion or mechanical defects must not be used for centrifuging.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Hettich company are allowed to be utilised.
- This centrifuge is classified in Germany as a Group 3 device according to the *Medizinische Geräteverordnung MedGV* (the regulations on medical equipment).
- It conforms to safety regulations based on: IEC 61010-1, IEC 61010-2-020
- The safe operation and reliability of the centrifuge can only be guaranteed if:
  - the centrifuge is operated in accordance with the operating instructions,
  - the electrical installation on the site where the centrifuge is installed conforms to the demands of IEC stipulations,
  - prescribed tests to BGV A1, BGR 261 are carried out by an expert.
- With centrifuges for robotic use please pay attention to the notes of the key operated switch.

**No claim under guarantee will be considered by the manufacturer unless the above instructions have been adhered to.**

#### 4 Explanations of symbols and references

##### 4.1 Symbol used on the machine

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

##### 4.2 Symbols used in the operating instructions

This symbol refers to important circumstances.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

#### 5 Delivery checklist

- 1 Connecting cable
- 2 Fuses
- 1 Hexagonal pin-type spanner
- 1 Operating instructions
- 1 Notes on moving the equipment safely

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

#### 6 Unpacking the centrifuge

- Lift the carton upward and remove the padding.
- Lift the centrifuge on both sides with an appropriate number of helpers and place it on the laboratory table.

Do not lift by the front panel.

## 7 Initial operation

- According to the laboratory instrument standards EN 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system.  
This switch has to be placed remote from the centrifuge, prefered outside of the room in which the centrifuge is installed or near by the exit of this room.
- Remove the transportation safety device from the bottom of the housing, see sheet "Transportation safety device"
- Position the centrifuge in a stable and level manner in a suitable place. When setting up the equipment, care should be taken to provide the required safety area of 300 mm around the centrifuge in accordance with IEC 61010-1.

No persons and / or dangerous materials may be located within the safety zone when the centrifuge is running.

- Do not place any object in front of the ventiduct.  
Keep a ventilation area of 300 mm around the ventiduct.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Appendix, Technical data/Technical specification".
- Turn on the mains switch. Switch position "I".  
The machine type and program version will be displayed and the LEDs light up. After 8 seconds  
**◀ OPEN ▶ OEFFNEN** will be displayed on centrifuges with cooling and the left LED in the key **STOP / OPEN** will flash. On centrifuges without cooling the cover will open automatically and the last used centrifuge data will be displayed.
- Open the cover on centrifuges with cooling.  
The last used centrifuge data will be displayed.
- Remove the transportation safety device in the centrifuge, see sheet "Transportation safety device".

## 8 Opening and closing the lid

### 8.1 Opening the lid

- Press the button **STOP / OPEN**. The lid opens automatically and the left LED in the button **STOP / OPEN** go off.

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".


### 8.2 Closing the lid

- Place the lid and lightly press down the front edge of the lid. The locking action is effected by motor. The left LED in the button **STOP / OPEN** lights up.

Do not put your fingers between lid and housing.  
Do not bang the lid shut.

## 9 Installation and removal of the rotor

- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The motor shaft dog (D) has to fit in the rotor slot (B). The alignment of the groove is labelled on the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.


Check the firm seating of the rotor on a weekly basis.

- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.

## 10 Loading the rotor

- The rotors and hangers may only be loaded symmetrically. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- Always fill the centrifuge containers outside of the centrifuge.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.
- The weight of the permissible filling quantity is specified on each rotor. This weight may not be exceeded.

## 11 Aerosol tight sealing of angle rotors

To ensure aerosol sealing the lid of an aerosol sealed angular rotor must be firmly closed.

This is performed with the help of the supplied spanner that is inserted through the borehole in the rotary handle. If the rotary handle does not have a borehole then the lid must be firmly closed by turning in clockwise direction.

For available aerosol tight angle rotors, please see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

## 12 Control and display elements

Please turn to the rear side of this cover sheet.

Fig. 2, Fig. 3: Display and control panel

### 12.1 Control knob

For setting the individual parameters.

Turning anticlockwise reduces the value. Turning clockwise increases the value.

### 12.2 Control panel pushbuttons (keys)

- Selection control key for selection of specific parameter.

The subsequent parameter is selected by every further keystroke.

- Start centrifugation run. The LED in the button lights up during the centrifugation run as long as the rotor is turning.

- Short-term centrifugation.

The centrifugation run is effected as long as the button is held down. The LED in the button lights up during the centrifugation run as long as the rotor is turning.

- Store inputs and changes.

- End centrifugation run.
- The rotor runs down with a pre-selected brake step. The right-hand LED in the button lights up until the rotor is stationary. Once the rotor is stationary the left-hand LED flashes in the button. Pressing the button twice triggers the EMERGENCY STOP.

- Open lid.

The left-hand LED in the button goes out.

- Leave the parameter input.

- Switch between RPM and RCF display.  
RCF values are displayed in > <.

- Start pre-cooling.

The pre-cooling speed is settable. It is pre-adjusted to 10000 RPM.

### 12.3 Adjustment possibilities

**PROG RCL** Program position of the called-up program.

t/min	Running time. Settable from 0 - 99 min, in 1 min increments.
t/sec	Running time. Settable from 0 - 59 s, in 1 second increments.
	Continuous run ( $\infty$ ). Set parameter <b>t/min</b> and <b>t/sec</b> to zero.
RPM	Revolutions per minute. A numerical value from 500 RPM up to the maximum speed of the rotor can be set. Maximum speed of the rotor, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". Settable in increments of 10.
RAD/mm	Centrifugation radius. Input in mm. For centrifugation radius see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". The input of the radius is only possible if the RCF display (> RCF <) is selected.
RCF	Relative Centrifugal Acceleration. A numerical value can be set, which gives a speed between 500 RPM and the maximum speed of the rotor. Adjustable up to 100 in intervals of 1, and from 100 in intervals of 10. The RCF value is automatically rounded up or rounded down with regard to the RPM interval. The input of the RCF is only possible if the RCF display (> RCF <) is selected.
	Starting steps 1 - 9. Step 9 = shortest starting time, Step 1 = longest starting time.
	Brake steps 1 - 9. Step 9 = shortest run-down time, Step 1 = longest run-down time.
T/°C	Temperature Set Point (only in centrifuges with cooling). Adjustable from -10°C to +40°C, in 1°C intervals. The lowest obtainable temperature depends on the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
<b>PROG STO</b>	Program position on which the program is stored. 4 programs can be stored (program positions 1 - 2 - 3 - 4).

## 13 Programming

### 13.1 Programme input/alteration

 If no key is pressed for 8 seconds long after the selection or during the input of parameters, the previous values will be shown in the display. The input of parameter then has to be executed again.

- Select the RPM or RCF display with the key **RCF**. RCF values are displayed in **> <**.
- Select the desired parameters using the button **SELECT** and set using the knob **○**. If continuous run is desired the parameters t/min and t/sec need to be set to zero on the knob **○** and the symbol  $\infty$  is displayed.
- The parameter **PROG STO** can be selected using the button **SELECT** and the desired program position set using the knob **○**.
- Press the button **START / IMPULS** in order to store the setting on the desired program position. **\*\*\* ok \*\*\*** is displayed briefly as confirmation.

 The previous data in the program position is overwritten during saving.

### 13.2 Programme recall

- Select the parameter **PROG RCL** using the button **SELECT** and set the desired program position using the knob **○**.
- Press the button **START / IMPULS**. The centrifugation data of the selected program position is displayed.
- The parameters can be checked by pressing the button **SELECT**. To leave the parameter display press the button **STOP / OPEN** or press no button for a period of 8 seconds.

## 14 Centrifugation

 If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, the unbalance display lights up, and **IMBALANCE** is displayed.  
A centrifugation run can be stopped at any time by pushing the key **STOP / OPEN**. All parameters can be selected and altered during the centrifugation run (see Chapter "Programming"). You can switch-over at any time between the RPM and RCF display with the key **RCF**. The input of the centrifugation radius is necessary if you are working with the RCF display.  
If **OPEN OEFFNEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

- Turn on the mains switch. Switch position **I**.
- Load the rotor and close the centrifuge lid.

#### 14.1 Centrifugation with pre-set time

- Adjusting time or recall a programme with pre-set time (see Chapter "Programming").
- Press the key **START / IMPULS**. The LED in the button **START / IMPULS** lights up for as long as the rotor turns.
- After expiration of the time or with truncation of the centrifugation run by pushing the key **STOP / OPEN**, the run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the remaining time will be displayed. After the centrifuge has ended its run and the rotor has come to a halt, the lid opens automatically on centrifuges without cooling.

#### 14.2 Continuous run

- Adjusting the symbol  $\infty$  or recall a continuous run programme (see Chapter "Programming").
- Press the key **START / IMPULS**. The LED in the button **START / IMPULS** lights up for as long as the rotor turns. The time metering begins at 00:00.
- Press the key **STOP / OPEN** in order to stop the centrifugation run. The run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the expired time will be displayed. After the centrifuge has ended its run and the rotor has come to a halt, the lid opens automatically on centrifuges without cooling.

#### 14.3 Short-term centrifugation

- Hold down the key **START / IMPULS**. The LED in the button **START / IMPULS** lights up for as long as the rotor turns. The time metering begins at 00:00.
- Let go of the key **START / IMPULS** again in order to stop the centrifugation run. The run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the expired time will be displayed. After the centrifuge has ended its run and the rotor has come to a halt, the lid opens automatically on centrifuges without cooling.

### 15 Emergency Stop

- Press the key **STOP / OPEN** twice.

With Emergency Stop the run-down is effected with brake step 9 (shortest run-down time). Brake step 9 is displayed.

### 16 Acoustic Signal

The acoustic signal sounds:

- Upon the appearance of a disturbance in 2 second intervals.
- After completion of a centrifugation run and rotor standstill in 30 second intervals.

The acoustic signal is stopped by opening the lid or pressing any key.

The signal after completion of the centrifugation run can be activated or deactivated in the following manner, if the rotor is at standstill:

- Hold down the key **SELECT** for 8 seconds.  
After 8 seconds, **SOUND / BELL** appears in the display.
- Set using the knob  $\odot$  **OFF** or **ON**.
- Press the key **START / IMPULS** in order to store the setting.  
As confirmation, **\*\*\* ok \*\*\*** will be displayed for a short period.

### 17 Recall hours of operation

Recall hours of operation is only possible during rotor standstill.

- Hold down the key **SELECT** for 8 seconds.  
After 8 seconds, **SOUND / BELL** appears in the display.
- Press the key **SELECT** once again.  
The centrifuge's hours of operation (**CONTROL:**) are displayed.
- Press the key **STOP / OPEN** to exit the hours of operation recall.

## 18 Cooling (only in centrifuges with cooling)

The temperature set-point can be adjusted from -10°C to +40°C. The lowest obtainable temperature is dependent on the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").

### 18.1 Standby-cooling

With rotor standstill and closed lid the centrifugal chamber is cooled to the pre-selected temperature. The temperature set-point is shown in the display.

Standby cooling will be subject to a timed delay after a centrifuge run and the display will show **◀ OPEN  
◀ OEFFNEN**. The delay time can be pre-set in 1 minute steps from 1 to 5 minutes. It is pre-set to 1 minute.

With the rotor standing still and the cover open the delay time can be set as follows:

- Hold down the key **☒** for 8 seconds.  
After 8 seconds, **t/min = X** appears in the display.
- Use the rotary button **○** to set the delay time.
- Press the key **START / IMPULS** in order to store the setting.  
As confirmation, **\*\*\* ok \*\*\*** will be displayed for a short period.

To leave the delay time display press the key **STOP / OPEN** or do not press any key for a period of 8 seconds.

### 18.2 Pre-cooling the rotor

- Press the key **☒**. The LED in the button **START / IMPULS** lights up for as long as the rotor turns.
- Press the button **STOP / OPEN** to end the pre-cooling. The run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature and the expired time will be displayed.

The pre-cooling speed can be adjusted in decadic steps from 500 RPM to the max RPM of the rotor. It is pre-adjusted to 10000 RPM.

When the rotor is stationary and the lid open the pre-cooling speed can be set in the following manner:

- Hold down the key **☒** for 8 seconds.  
After 8 seconds, **t/min = X** appears in the display.
- Press the key **☒** once again.  
The set pre-cooling RPM - **RPM = XXXX** will be displayed.
- Set the desired pre-cooling speed using the knob **○**.
- Press the key **START / IMPULS** in order to store the setting.  
As confirmation, **\*\*\* ok \*\*\*** will be displayed for a short period.

To leave the pre-cooling RPM display press the key **STOP / OPEN** or do not press any key for a period of 8 seconds.

## 19 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left( \frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

**RCF** = relative centrifugal force

**RPM** = rotational speed (revolutions per minute)

**r** = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter " Anhang/Appendix,  
Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

## 20 Centrifugation of materials with higher density

The rotors are designed to centrifuge substances up to a maximum mean homogenous density of 1.2 kg/dm<sup>3</sup> when rotating at the stated speed.

Denser substances must be centrifuged at lower speed.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}\text{)} = \sqrt{\frac{1.2}{\text{Greater density}}} \times \text{Rated speed}$$

e.g.: RPM 4000, density 1.6 kg/dm<sup>3</sup>

$$\text{n}_{\text{red}} = \sqrt{\frac{1.2}{1.6}} \times 4000 = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

## 21 Emergency release

In the event of a power failure the lid cannot be opened by motor. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.  
Open the lid only during rotor standstill.

Please turn to the rear side of the cover sheet.

- Insert the hexagonal wrench key into the bore hole (Fig. 1, A) and carefully rotate by half a turn in clockwise direction until the lid can be opened.
- Pull the hexagon socket head wrench out of the drilling again.
- After turning the centrifuge on again, press the button so that the motor-driven lid locking once again assumes the normal position (opened).

## 22 Maintenance and servicing

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Cleaning agents and disinfectants which lie in the pH range 5 – 8 are to be utilised. Alkaline cleaning agents with a pH value > 8 are to be avoided.
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

### 22.1 Centrifuge

- Regularly clean the centrifuge housing and the centrifugal chamber; and if necessary, clean with soap or a mild cleaning agent and water. For one thing, this serves purposes of hygiene, and it also prevents corrosion through adhering impurities.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately. For surface disinfection we recommend Bacillol Plus from the company Bode Chemicals in Hamburg. After use the disinfection medium Bacillol Plus must be thoroughly removed with water.
- Lightly grease the rubber seal of the centrifugal chamber after every cleaning.

## 22.2 Rotors and Attachments

- Rotors and accessory parts must be regularly cleaned with soap or a mild cleaning agent and water in order to prevent corrosion and changes of material. Cleaning is recommended at least once a week, even better after every usage.
- If the rotor or accessory parts are contaminated by pathogenic or radioactive material, a suitable cleaning has to be executed. For disinfection we recommend Heliipur H plus N from the company B. Braun Melsungen. For the removal of radioactive material we recommend deco neutracon from the company Decon Laboratories Limited.
- The rotors and accessory parts must be dried immediately after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- With aerosol-sealed rotors and bio safety systems (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories") the sealing rings are to be checked and cleaned regularly (weekly). The sealing ring is to be replaced immediately upon indication of crack formation, embrittlement or abrasive wear.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and accessory parts are to be checked on a monthly basis for corrosion damage.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

## 22.3 Autoclaving

Swing-out rotors, angle rotors made of aluminium and accessory parts can be autoclaved at 121°C (20 min).

The lids of the rotors and containers must be removed prior to autoclaving.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

## 22.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splinters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

## 23 Faults



If the fault cannot be rectified according to the faults table, Hettich customer services must be informed.

Please state the type of centrifuge and the factory serial number. Both values are visible on the centrifuge type plate.

MAINS-RESET: - Mains switch OFF, for longer than 10 seconds.  
- Mains switch ON.

Message / fault		Cause	Remedy
No display	---	- No voltage. - Mains input fuses defective.	- Check supply voltage. - Check mains input fuses. - Mains switch ON.
TACHO - ERROR	01	- Faulty speedometer.	- Open lid.
	02	- Motor, electronics defective.	- Turn rotor manually. - MAINS-RESET when power is switched on, rotor should turn.
IMBALANCE	---	Imbalance about motor axis through weight differential in rotor assembly.	- Open lid. - Correct imbalance.
CONTROL - ERROR	04	Error in lid locking or lid closure.	- Open lid.
N > MAX	05	Rotation too fast	- MAINS-RESET
N < MIN	13	Rotation too slow	
MAINS INTERRUPT	---	Power failure, centrifuging not properly completed	- Open lid. - Push  button.
VERSIONS-ERROR	12	Mismatch between electronic components	- MAINS-RESET
CONTROL-ERROR	21, 22, 24 - 27	Error / defect electronics	
CONTROL-ERROR	23	Error / defect controls	
SER I/O - ERROR	30 - 38	Error / defect electronics	
° C * - ERROR	50 - 56	Error / defect electronics	
FU / CCI - ERROR	60 - 83	Error / defect electronics / motor	
SYNC-ERROR	90	Error / defect electronics	
SENSOR-ERROR	91, 92	Error / defect unbalance sensor	
KEYBOARD-ERROR	---	Error / defect controls	

## 24 Change mains input fuse

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuse.
- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

## 25 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

You will find further information on our homepage: <http://www.HettichLab.com>

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All rights reserved. No part of this publication may be reproduced without the written prior permission of the copyright owner.