
Installation Instructions | Betriebsanleitung | Mode d'emploi | Istruzioni per l'uso |  
Instrucciones de instalación | Handleiding | Instruções de funcionamento | Bruksanvisning |  
Käyttöohje | Betjeningsvejledning | Εγχειρίδιο οδηγιών χρήσης

## Sartorius Combics

CAPXS.. Models

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

98648-010-96

# Contents

**English page 2**  
In cases involving questions of interpretation, the German-language version shall prevail.

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- 2 Symbols
- 2 Intended Use
- 3 Warning and Safety Information
- 3 Warranty
- 4 Installation Instructions
- 6 Installation
- 9 Care and Maintenance

## The following symbols are used in these instructions:

- indicates required steps
- indicates steps required only under certain conditions
- > describes what happens after you have performed a certain step
- indicates an item in a list
- △ indicates a hazard

Make sure you observe the following warning and safety information in its entirety during installation and operation, as well as while performing maintenance and repair work on the equipment. It is important that all personnel using the Combics equipment understand this information, and have access to it at all times.  
Furthermore, the warning and safety information supplied with any electrical equipment connected, such as the indicator, must be observed as well. The warning and safety information can be supplemented by the equipment operator.  
Make sure all operating personnel are informed of any additions to these instructions.

## Intended Use

The weighing platform and the connected indicator are intended exclusively for use in weighing.

# Warning and Safety Information

The weighing platform meets the requirements for Group II, Category 2 equipment in accordance with EC Directive 94/9/EC and bears the designation  $\text{Ex II 2 GD EEx ib IIC T4...T6}$  T135...155°C in accordance with KEMA EC type-examination certificate 02ATEX1010X.

Furthermore, the weighing platform meets the EC Directives for electromagnetic compatibility (see the Declaration of Conformity). Improper use or handling, however, can result in damage and/or injury.

- The weighing platform can be operated indoors or outdoors. This weighing platform may not be used in the medical industry (as a medicinal product).
- The weighing platform may be used in Zone 1, 2, 21 or 22 hazardous areas or Class 1, Division 1 or 2 hazardous locations (gas and dust explosion hazards). Please make sure the currently valid regulations and guidelines for installing equipment in the hazardous areas/locations listed above (e.g., EN60079-14 or, in Germany, ExeV) are strictly observed. Whether the equipment can be used in a given area containing potentially explosive agents must be checked on a case-by-case basis.
- The weighing platform may be used and operated by qualified personnel only. The permitted uses of the weighing platform are specified in the type-examination certificate.
- Do not expose the weighing platform to aggressive chemical vapors or to extreme temperatures, moisture, shocks, or vibration. The permissible operating temperature range during operation is -10°C to 40°C (14°F to 104°F).
- If you use suction lifting equipment to lift the load plate, always wear gloves, hard-toed safety boots and protective clothing. Warning: Danger of personal injury! Only reliable personnel who are qualified to perform such work are allowed to use suction lifting equipment.
- Suspension points are designated on models of 1000 x 1000 mm and larger. If you need to transport or lift the scale or load plate using a crane, do not stand underneath the suspended scale or load plate.
- Make sure to observe the applicable safety rules and regulations for the prevention of accidents.
- Do not damage the junction box or the load cells during transport.
- Installation in a Zone 1, 2, 21 or 22 hazardous area or Class 1, Division 1 or 2 hazardous location must be performed by a trained technician who is familiar with the assembly and operation of the equipment, as well as with

the procedure for putting the system into operation. Furthermore, the trained technician must have the required qualifications and must be familiar with the relevant guidelines and regulations. If you need assistance, contact your Sartorius dealer or the Sartorius Service Center. Any installation work that does not conform to the instructions in this manual will result in forfeiture of all claims under the manufacturer's warranty. Be sure to observe all restrictions listed in the type-examination certificate. Operating the weighing platform beyond the limits imposed by these restrictions is not permitted, and constitutes use of the equipment for other than its intended purpose.


- Have the equipment inspected at appropriate intervals for correct functioning and safety by a trained technician.
- The junction box may be opened only by authorized service technicians who have been trained by Sartorius and who follow Sartorius' standard operating procedures for maintenance and repair work.
- Always make sure the weighing platform is disconnected from AC power before performing any installation, cleaning, maintenance or repair work. If the equipment housing is opened by anyone other than persons authorized by Sartorius, all claims under the manufacturer's warranty are forfeited. Use only original Sartorius spare parts.
- Handle the equipment as specified in EN 60529 (in Germany: VDE 0470 Part 1) in accordance with its protection (IP) rating. The protection rating of the weighing platform is IP67. Do not damage the IP protection when cleaning the equipment. The IP protection rating is ensured only if the rubber gasket is installed on the junction box and all cable gland screw fasteners are connected securely. Any installation work that does not conform to the instructions in this manual will result in forfeiture of all claims under the manufacturer's warranty.
- If you use cables purchased from another manufacturer, check the pin assignments before connecting the cable to Sartorius equipment. Check the assignments in the cable against those specified by Sartorius and disconnect any wires that are assigned differently. The operator shall be solely responsible for any damage or injuries that occur when using cables not supplied by Sartorius.

- When using the weighing platform in hazardous areas/locations, make sure there is no current or voltage in the equipment before connecting or disconnecting current-carrying cables to or from the platform. Disconnect the platform from AC power before connecting or disconnecting cables.
- Avoid exposing the weighing platform to static electricity; be sure to connect the equipotential bonding conductor to the junction box.
- Disconnecting equipotential bonding conductors is not permitted.
- If you see any indication that the weighing platform cannot be operated safely (for example, due to damage), turn off the platform and lock it in a secure place so that it cannot be used for the time being. Observe the relevant safety precautions and inform personnel as required.
- The casing on all connecting cables, as well as the casing on wires inside the equipment housing, is made of PVC. Chemicals that corrode these materials must be kept away from these cables.
- Make sure the weighing instrument is not exposed to substances that release chlorine ions at the place of use. If such exposure cannot be ruled out, the operator is responsible for establishing and observing appropriate safety precautions, to be checked at regular intervals for continued effectiveness.
- The load cells installed are manufactured by Global Weighing Technologies; type designations: 011299/... through 011311/...; EC type-examination certificate no. PTB 02 ATEX 2061. The designation code is:  $\text{Ex II 2 G EEx ia IIC T4/T6}$ . The specifications are listed in the EC type-examination certificate. The load cell may be used and operated by qualified personnel only. The permissible use of the equipment is defined by the equipment specifications and the relevant safety regulations. Operating the load cell beyond the specifications listed in the type-approval certificate is not permitted, and constitutes use of the equipment for other than its intended purpose. All restrictions listed in the type-examination certificate must be strictly observed.


## Warranty

Do not miss out on the benefits of our full warranty. Please contact your local Sartorius office or dealer for further information. If available, complete the warranty registration card, indicating the date of installation, and return the card to your Sartorius office or dealer.


# Installation Instructions


- Prepare a suitable place of installation for the weighing platform. The place of installation should be dry, level and even. The allowable operating temperature range is  $-10^{\circ}\text{C}$  to  $+40^{\circ}$  ( $14^{\circ}\text{F}$  to  $104^{\circ}\text{F}$ ). The permissible load that can be carried by the chosen working surface must be sufficient for both the weighing platform and any load placed on the platform.  
If you need to use the weighing platform in areas exposed to heavy traffic (e.g., forklift trucks), you should install a protective frame, consisting of angular braces, around the weighing platform.  
Do not expose the weighing platform unnecessarily to extreme temperatures, moisture, shocks, or vibration that could result in damage.


- The air bubble must be centered within the circle on the level indicator.

- Level the weighing platform using the leveling feet as described below:
  - Check to ensure that all leveling feet rest securely on the work surface.
    - > Each of the leveling feet must support an equal load.
  - Loosen the lock nuts on the leveling feet using a 19-mm open-end wrench (spanner).
    - > Adjusting the leveling feet:  
To raise the weighing platform, extend the leveling feet (turn clockwise).  
To lower the weighing platform, retract the leveling feet (turn counterclockwise).
  - After leveling the weighing platform, retighten the lock nuts securely as described below:  
Low-capacity platforms (1 load cell): tighten the lockouts against the platform frame;  
high-capacity platforms (4 load cells): tighten the lock nuts against the platform feet.

- If the weighing platform is in a hazardous area/location, it must be grounded (i.e., an equipotential bonding conductor must be connected). This connection should be made by a trained technician.  
All Combics weighing platforms are equipped with a connector for the grounding conductor.  
This is located either below the load pan, on the junction box, or on the lower frame of the weighing platform. The position is marked in each case by the symbol shown here, indicating the grounding (earthing) connection.  
The grounding conductor is connected to a threaded bolt or terminal screw, or a bore hole is provided. If a bore hole is provided, use a stainless steel screw and nut to connect the grounding conductor. Use of a tooth lock washer is recommended, to prevent the screw from coming loose. The wire used for the grounding conductor should have a cross-sectional diameter of at least  $4\text{ mm}^2$ , with a suitable ring lug attached.  
Connect all equipment, including peripheral devices, to the equipotential bonding conductor.

### Conditions for Installation

Before putting the equipment into operation, it is important to make sure that the cable of the power supply is correctly connected to the power outlet (mains supply). All equipment must be connected to the equipotential bonding conductor via grounding cable (not included in delivery) connected to the grounding terminals on each device. The dimensions of the grounding cable are specified in national regulations for electrical installations. Installation must be performed by a trained technician in accordance with national regulations and generally acknowledged rules of engineering.

Use only cabling and extensions approved by Sartorius, as these are made in accordance with the restrictions on permissible cable lengths imposed by both the capacitance and inductivity values (see the "Annex to EC Type-Examination Certificate" in the operating instructions) and the requirements for electromagnetic compatibility.

Before putting the weighing system into operation for the first time, make sure there is no hazard of explosion present at the place of installation. If there is any indication that the equipment does not function properly (e.g., display remains blank, or no display backlighting) due to damage during transport, disconnect the equipment from power and notify your nearest Sartorius Service Center.


The weighing platform specifications for Ui, Li, Pi, Ta, temperature class, Ci, and Li are listed in the EC type-examination certificate. This certificate also specifies the indicator models that can be used with these platforms. These specifications must be observed when connecting an indicator to the platform. Ci and Li values apply to weighing platforms with a 3-m cable for connecting the indicator. (See also the "Certificate of Proof of Intrinsic Safety" in the appendix to this manual.)

The explosion-protected weighing system must be installed in accordance with acknowledged rules of engineering. These include national and international laws and regulations, such as EN 60079-14 and EN 5281-1-2, for example. In particular, the conditions described under Item 17 of the KEMA EC type-examination certificate, "Special Conditions for Safe Use", must be observed. Furthermore, national regulations for accident prevention and environmental protection must be observed at all times.

Before the platform is operated in a hazardous area/location, it must be inspected either by a certified electrician or under the guidance and supervision of a certified electrician to make sure that the weighing system complies with the applicable regulations (in Germany, Section 12 of the ExelxV). Determine whether the installation must be registered with technical inspection authorities (such as the trade board) in your country. Regular inspections must also be performed on the system during operation. The system should be inspected at intervals short enough to permit the prevention or early detection of defects that arise as a result of normal wear and tear. The longest permissible interval period is 3 years. Other conditions and standards that regulate the installation and operation of the equipment and are applicable in your country must be met as well. When performing inspections, generally acknowledged rules of engineering relevant to these conditions must also be applied.

If the terminal housing is opened by anyone other than persons authorized by Sartorius, or if the terminal is installed or operated incorrectly, this will result in forfeiture of the approval for use in the stated hazardous area(s)/location(s) and of all claims under the manufacturer's warranty.


# Installation

- Connect the cable of the weighing platform to a suitable indicator, such as the FCT01-X from Sartorius.

Note:

The cable gland on Sartorius indicators is installed at the factory. Please use extreme caution when performing any work on the equipment that affects this cable gland.


- Remove the casing from a length of the cable end and install the cable as follows:
  - Remove the screw fastener from the cable gland and slide it over the end of the cable.
  - Remove the casing from a section of the cable end (see illustration). The shielding (1) must have contact with the clamps (2). Fold back 3 to 4 mm of the shielding so that it overlaps the O-ring (see illustration).
  - Expose approximately 15 cm (4 inches) of the wires (3) for connection within the terminal.
  - Thread the cable through the cable gland.
  - Make sure the shielding is in contact with the clamps. The shielding supplies the grounding connection.
  - Replace the screw fastener and tighten it securely.

- Connect the wire ends to the terminal strip:

- Remove the casing from a section of the cable end. Expose enough of the wires in the cable for installation.
- Remove the casing from approx. 1 cm (0.5 in.) of the wires and affix ferrules to the wire ends.
- Slide the enclosed protective sleeve over all the wires down to the cable gland, so that no more than 5 cm of the wires is exposed between the end of the sleeve and the terminal strip.
- Attach the wires securely to the screw terminals (blue = positive, brown or black = negative)

## Color Codes in the Weighing Platform Connecting Cable, Model CAPXS..

No.	Signal des.	Meaning	Color of wire casing
1	BR_POS	Bridge supply voltage (+)	blue
2	SENSE_POS	Sense (+)	green
3	OUT_POS	Bridge supply voltage	white
4	OUT_NEG	Measuring voltage positive	red
5	SENSE_NEG	Measuring voltage (-)	gray
6	BR_NEG	Sense (-)	black or brown

- Use a screwdriver to tighten the terminal screws.

## Note on connecting the FCT01-X:

The bridge supply voltages are protected by 62mA fuses. To prevent the fuses from reacting if wires are incorrectly connected, resistors between pin 1 and pin 6, between pin 1 and the housing and between pin 6 and the housing must measure greater than 140 ohms. Low-ohm connections are also required between pins 1 and 2 and between pins 5 and 6.

**Key to Model Codes (CAPX ab - cd - e):**

Product family	Material	Number of load cells	Capacity in kg see Table 2	Dimensions see Table 1	Resolution see Table 2
CAPX..	a	b	c	d	e
Combics	S =	1 = one load cell			
Analog Platform	stainless steel	4 = four load cells	3 6 15 30 60 150 300 600 1500 3000		

Table 1: Model-specific data: Dimensions

Designation	DC	ED	FE	GF	IG	LL	NL	RN	RR	WR
Width (mm)	240	300	400	500	600	1000	1000	1250	1500	1500
Length (mm)	320	400	500	650	800	1000	1250	1500	1500	2000


**Cable lengths (power cord)**

Length (m), approx.	1.5	1.5	3.0	3.0	3.0	6.0	6.0	6.0	6.0	6.0

Table 2: Resolutions:

1 weighing range			2 weighing ranges*						
	-L	-I	-LCE		-NCE	2 × 3000e			
Weighing range	15,000d	30,000d	1×3000e		Weighing range	Weighing range 1	Resolution range 1	Weighing range 2	Resolution range 2
in kg	in g	in g	in g		in kg	in kg	in g	in kg	in g
3	0.2	0.1	1		3	1.5	0.5	3	1
6	0.5	0.2	2		6	3	1	6	2
15	1	0.5	5		15	6	2	15	5
30	2	1	10		30	15	5	30	10
60	5	2	20		60	30	10	60	20
150	10	5	50		150	60	20	150	50
300	20	10	100		300	150	50	300	100
600	50	20	200		600	300	100	600	200
1500	100	50	500		1500	600	200	1500	500
3000	200	100	1000		3000	1500	500	3000	1000

- \* The weighing ranges allowed in legal metrology in the EU are listed in the Declaration of Conformity. On weighing instruments with two weighing ranges (2 × 3000e), the instrument does not automatically switch back to the higher resolution of the fine range once you use the second, higher-capacity range. In other words, the lower resolution of the higher range is retained.  
Example: CAPXS4-3000WR-NCE = Combics analog platform in stainless steel with 4 load cells, weighing capacity 3000 kg, 2 × 3000e in verifiable version.

#### Operating Limits

Do not place loads on the weighing platform that exceed its maximum weighing capacity.

Depending on the position of the load (center, side, one-sided corner load), the maximum capacity of the weighing platform is as follows:

Model	Center	Side	Corner
DC	50	35	20
ED	130	85	45
FE	300	200	100
GF	450	300	150

IG	900	600	300
LL	4500	3000	1500
NL	4500	3000	1500
RN	4500	3000	1500
RR	4500	3000	1500
WR	4500	3000	1500

⚠ Please follow the safety instructions!

# Care and Maintenance

## Cleaning

- Unplug the weighing instrument from the AC power before cleaning.  
Remove dust from the weighing platform regularly (see EC type-examination certificate)
- If the weighing instrument is in a dry room, use a damp cloth to wipe down the weighing platform. You can use common household cleaning agents. Follow the manufacturer's instructions for the cleaning agent.
- The leveling feet are covered with rubber buffers. In hazardous areas/locations, these must be cleaned using only a damp cloth. This prevents build-up of static electricity.
- △ Never use concentrated acids, bases, solvents or pure alcohol to clean the weighing platform.
- To clean the weighing platform in a wet area, wash it down using a gentle stream of water (60°C max.) sprayed over the top of the load plate.
- △ Do not use high-pressure cleaning equipment to clean the weighing platform.
- > If the water that you use to clean the weighing platform is too hot or too cold, the difference in temperature between the water and the weighing platform can cause condensation within the weighing platform. This condensation may cause the weighing platform to malfunction.
- If the scale is installed in a pit, make sure that no debris builds up in the crevices between the pit and the platform, to prevent weighing errors.
- Regularly remove debris from the bottom of the pit.

## Cleaning the Inside of the Weighing Platform

- To clean the inside of the weighing platform, remove the load plate.  
Be especially careful when removing the load plate from scales of 1000 x 1000 mm or larger.
- △ Please follow the safety instructions.
- Use compressed air to blow debris out of the inside of the scale or flush out using a gentle stream of water (60°C max.). Make sure that no debris builds up in the gap between the load receptor and the fastening plate, to avoid compromising the overload protection.

## Corrosive Environment

- Remove all traces of corrosive substances from the weighing platform on a regular basis.

## Instructions for Recycling the Packaging

If you no longer need the packaging after successful installation of the equipment, you should return it for recycling, as it is a valuable source of secondary raw material. Sartorius products are packaged to ensure safe shipment using environmentally friendly materials. For information on recycling options, including recycling of old weighing equipment, contact your municipal waste disposal center or local recycling depot. Make sure you remove the batteries before disposing of the equipment. Rechargeable batteries are toxic waste and as such may not be disposed of with normal household waste. Contact your municipal waste disposal center or local recycling depot for details on the proper disposal of batteries.

## Model-specific Information: Serial Number Coding

The month and year of manufacture are encoded in the serial number as follows:

Y M M x x x x

Y=	Year:
1	2000–2006
2	2007–2013
3	2014–2020
4	2021–2027
5	2028–2034
6	2035–2041
7	2042–2048
8	2049–2055
9	2056–2062

The first digit represents a 7-year period as indicated in the table above. The next 2 digits represent the month. The months are numbered consecutively, starting with 13, over the entire 7-year period. Thus the number representing the month also indicates the specific year of manufacture.

2000 13–24

2001 25–26

...

Example:

113xxxxx -> January 2000

The individual devices are numbered consecutively in the last 5 digits, starting from 00000 again at the beginning of each month.

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Copyright by Sartorius AG,  
Goettingen, Germany.  
All rights reserved. No part  
of this publication may  
be reprinted or translated in  
any form or by any means  
without the prior written  
permission of Sartorius AG.  
The status of the information,  
specifications and illustrations  
in this manual is indicated  
by the date given below.  
Sartorius AG reserves the  
right to make changes to  
the technology, features,  
specifications and design of  
the equipment without notice.

Status:  
April 2004, Sartorius AG,  
Goettingen, Germany

Distributore Autorizzato :  
Geass S.r.l. - Torino  
Tel.: +39 011.22.91.578  
[info@geass.com](mailto:info@geass.com)  
web site : [www.geass.com](http://www.geass.com)

Printed in Germany on paper that has been bleached  
without any use of chlorine  
MJ · W3A000 · KT  
Publication No.: WCA6003-u04042