
**Installation Instructions and Safety Information | Installation und Sicherheitshinweise**  
**Installation et conseils de sécurité | Istruzioni per l'installazione e di sicurezza**  
**Advertencias de Instalación y Seguridad**

## Sartorius EC1XS.-.....-L

Electronic Stainless Steel Scales for Use in Hazardous Areas  
Elektronische Edelstahlwaagen für explosionsgefährdete Bereiche  
Balances électroniques en acier inoxydable pour atmosphères explosives  
Bilance elettroniche in acciaio inox per le aree a rischio di esplosione  
Básculas electrónicas de acero inoxidable para las zonas con riesgo de explosión

98648-013-65

---

**English – page 3**

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
# General View of the Equipment


- 1 Display and control unit (evaluation unit)
- 2 Display
- 3 Print key
- 4 Function key
- 5 Tare key
- 6 Zero key
- 7 On/off/standby key
- 8 Weighing platform
- 9 Leveling foot
- 10 Level indicator

- 11 Data output port: optional (RS-422)
- 12 Ground/earth terminal
- 13 Data output port: standard equipment (RS-232)
- 14 Menu access switch
- 15 Connecting cable (weighing platform to display and control unit)
- 16 Connector for foot switch or other optional peripheral device
- 17 DC connector for power supply: YPS05-ZD (non-hazardous area) or YPS05-XD (hazardous area)


## Configuration Example

EC1XS.-.....-L complete scale with power supply and an interface converter for data transfer

Weighing platform


YPS05-XD display and control unit

YPS05-ZDR

YDI05-Z  
(interface converter)

Computer

Hazardous area  
Installation in zone 1/2/21/22

Non-hazardous area

## Contents

3	<b>General View of the Equipment</b>
4	<b>Intended Use</b>
4	<b>Warnings and Safety Precautions</b>
6	<b>Installation Instructions</b>
8	<b>Installation</b>
8	<b>Connecting the EC1XS.-.....-L</b>
9	<b>Care and Maintenance</b>
10	<b>Declarations of Conformity</b>
11	<b>Documents for Hazardous Areas</b>

### The following symbols are used in these instructions:

- indicates required steps
- indicates steps required only under certain conditions
- > describes what happens after you have performed a particular step
- indicates an item in a list
- △ indicates a hazard

Improper use or handling can result in damage and/or injury. The scale may be installed and operated by qualified personnel only. Make sure you observe the warning and safety information in its entirety during installation and operation, as well as while performing maintenance and repair work on the equipment. The regulations, applicable laws, occupational safety requirements and environmental protection laws valid in your country must be observed. In particular, the conditions described under Item 17 of the KEMA EC type-examination certificate, "Special Conditions for Safe Use," must be observed. It is important that all personnel using the equipment understand this warning and safety information, and have access to the relevant documents at all times. Furthermore, the warning and safety information supplied with any electrical equipment connected, such as peripheral devices, must be observed as well. The warnings and safety precautions may have to be supplemented by the equipment operator. All operating personnel must be informed of any additions to these instructions. Make sure the equipment is accessible at all times.

## Intended Use

The compact scales in the EC1XS-.....-L series are designed for mass determination in zone 1, 2, 21 and 22 hazardous areas. The scale consists of a weighing platform and a display and control unit (evaluation unit). An external power supply is required for operation of the scale (option YPS05-ZD for installation in the non-hazardous area, or YPS05-XD for the hazardous area). In addition to the basic weighing function, the scale is equipped with programs for toggling weight units and gross/net weighing. Weight values can be transferred through a Zener barrier (e.g., option YDI05-Z) over an RS-232 port (built-in) or an RS-422 port (optional) to a printer or computer outside the hazardous area. A foot switch can also be connected to the scale.

### The display and control unit (evaluation unit) comes with a separate operating manual.

#### Make sure all safety instructions are observed.

##### Note:

Read the installation and operating instructions carefully before connecting and operating the EC1XS.-.....-L. Keep the instruction manuals in a safe place.

## Warnings and Safety Precautions

### General Provisions for Installing Explosion-protected EC1XS.-.....-L Models in Zone 1, 2, 21 and 22 Hazardous Areas

Model EC1XS.-.....-L scales meet the requirements of EC Directive 94/9/EC for Group II, Category 2G or 2D devices and are labeled in accordance with EC type-examination certificate KEMA04ATEX1279 (see the scale operating instructions).

Furthermore, EC1XS.-.....-L models meet the requirements of EC Directives for electromagnetic compatibility and electrical safety (see "Declaration of Conformity" in the scale operating instructions).


- The area of use for EC1XS.-.....-L models and accessories (power supplies, interface converter) is defined in the type-examination certificate. All restrictions listed in the type-examination certificate must be observed. Operation outside of these restrictions is not permitted, and is considered use of the EC1XS.-.....-L for other than its intended purpose. Any installation work that does not conform to the instructions in this manual will result in forfeiture of all claims under the manufacturer's warranty.
- When the equipment is operated in a hazardous area, the national laws and regulations applicable at the place of use must be observed (e.g., EN60079-14 (gases); EN50281-1-2 (dusts)). Ask your supplier for information on the legal regulations applicable in your country.
- The EC1XS.-.....-L is may not be used in the medical industry (as a medicinal product).
- Whether the equipment can be used in a given area containing potentially explosive agents must be checked on a case-by-case basis.
- If the equipment housing is opened by anyone other than persons authorized by Sartorius, this will negate its conformity with regulations governing use in the stated hazardous areas and result in forfeiture of all claims under the manufacturer's warranty.
- Installation of the system must be performed by a certified electrician who is familiar with the assembly, start-up and operation of both the system and the relevant guidelines and regulations, and has the required qualifications for performing the installation. If you need assistance, contact your Sartorius dealer or the Sartorius Service Center.
- Avoid generating static electricity. The weighing platform is connected to the display and control unit via the shield in the connecting cable, or via an optional column.

- Connect all equipment and accessories to the equipotential bonding conductor; see also "Documents for Hazardous Areas." Each grounding terminal is marked by a "ground" symbol . Disconnecting equipotential bonding conductors is not permitted. Make sure the connection cannot come loose. The grounding conductor is connected to a threaded bolt or terminal screw, or a bore hole is provided. If a bore hole is provided, use a stainless steel screw and nut to connect the grounding conductor. Use of a tooth lock washer is recommended to prevent the screw from coming loose. The wire used for the grounding conductor should have a cross-sectional diameter of at least 4 mm<sup>2</sup> and have a suitable ring lug attached. If the weighing platform is moved on rollers in a hazardous area, the speed of movement must be limited to ≤1 m/s or the scale must be grounded.
  - Responsibility for ensuring correct and effective grounding of the equipment lies with the system operator.
  - When using the model EC1XS.-.....-L scale in a hazardous area, make sure there is no current or voltage in the equipment before connecting or disconnecting current-carrying cables to or from the scale or system. Disconnect the system from AC power before connecting or disconnecting cables.
- The specifications (see EC type-examination certificate) must be taken into account when connecting accessory equipment to the platform. See also "Documents for Hazardous Areas."
- Any sockets or ports on the display and control unit to which no cable is connected must be protected within hazardous areas by affixing a cap with the required IP protection rating.
  - Before the EC1XS.-.....-L scale is operated in a hazardous area, it must be inspected either by a certified electrician or under the guidance and supervision of a certified electrician to make sure that the system complies with the applicable regulations. Determine whether the installation must be registered with the technical inspection authorities (such as the trade board) in your country. Regular inspections must be performed on the system during operation. The system should be inspected at intervals scheduled to allow for the prevention or early detection of defects that arise as a result of normal wear and tear. The longest permissible interval period is 3 years. Other conditions and standards that regulate the safety of installation and operation of the equipment and are applicable in your country must be met as well. When performing inspections, generally acknowledged rules of engineering relevant to these conditions must be applied.
  - Do not expose the scale to aggressive chemical vapors or to extreme temperatures, moisture, shocks, or vibration. The allowable ambient temperature range during operation is -10°C to +40°C (14°F to +104°F). Make sure the place of installation is adequately ventilated to prevent build-up of excessive heat.
  - If you use cables purchased from another manufacturer, check the pin assignments in the cable against those specified by Sartorius before connecting the cable to Sartorius equipment, and disconnect any wires that are assigned differently. The operator shall be solely responsible for any damage or injuries that occur when using cables not supplied by Sartorius. Use only genuine Sartorius spare parts.
- #### For the User
- Always make sure the equipment is disconnected from AC power before performing any installation, cleaning, maintenance or repair work on the scale or weighing system.
  - If you see any indication that the system cannot be operated safely (for example, due to damage), turn it off and lock it in a secure place or otherwise prevent use of the equipment for the time being.
  - Chemicals (such as gases or dusts) that can corrode and damage the inside or outside of the device must be kept away from the equipment. Handle the equipment and accessories in accordance with the protection (IP) rating (EN 60529).
  - The casing on all connecting cables, as well as the casing on wires inside the equipment housing, is made of PVC. The casing on the power cable is made of rubber.

# Installation Instructions

## Unpacking the Equipment

- After unpacking the equipment, please check it immediately for any visible damage. Please read the enclosed operating instructions carefully before putting the equipment into operation.
- If damage is evident, refer to the instructions under "Safety Inspection" in the chapter entitled "Care and Maintenance."
- It is a good idea to save the box and all parts of the packaging until you have successfully installed your equipment. Only the original packaging provides the best protection for shipment.
- Before packing your equipment, unplug all connected cables to prevent damage.


## Equipment Supplied

- Complete scale (weighing platform and display and control unit)
- Installation Instructions and Safety Information (incl. Documents for Hazardous Areas)
- Operating instructions (incl. approvals for optional equipment)  
An external power supply is necessary for operation of the scale:  
YPS05-XD (hazardous area)  
YPS05-ZD (non-hazardous area)


### Note:

No power supply is included with the equipment.

## Requirements on the Place of Installation

- Observe all safety instructions. The EC1XS.-.....-L is designed to provide reliable results under normal ambient conditions.  
When choosing a location to set up the EC1XS.-.....-LX, observe the following so that you will be able to work with added speed and accuracy:
  - Set up the EC1XS.-.....-L on a stable, even surface. The carrying capacity of the chosen working surface must be sufficient for both the weighing platform and any load placed on it.
  - Avoid placing the equipment in close proximity to a heater or otherwise exposing it to heat or direct sunlight.
  - protect the EC1XS.-.....-L from direct exposure to drafts that come from open windows or doors.
  - Avoid exposing the EC1XS.-.....-L to excessive vibration.

- Protect the EC1XS.-.....-L from aggressive chemical vapors.
- Do not expose the equipment to excessive electromagnetic interference, as it can affect the value displayed.
- Turn off the power when the system is not in use
- The allowable operating temperature range is -10°C to +40°C (14°F to 104°F)


- Level the weighing platform using the leveling feet as described below:
- Check to ensure that all leveling feet rest securely on the work surface.
  - > Each of the leveling feet must support an equal load.
- Loosen the locknuts on the leveling feet using a 19-mm open-end wrench (spanner).
- > Adjusting the leveling feet:
  - To raise the weighing platform, extend the leveling feet (turn clockwise).
  - To lower the weighing platform, retract the leveling feet (turn counterclockwise).
- After leveling the weighing platform, tighten the lock nuts as follows:  
Low-capacity platforms (1 load cell): tighten the locknuts against the platform frame. High-capacity platforms (4 load cells): tighten the locknuts against the platform feet.
- If the scale is in a hazardous area, it must be grounded (i.e., an equipotential bonding conductor must be connected; resistance: ≤1 ohm). This connection should be made by a trained technician (see "Warnings and Safety Precautions").

⚠ Sartorius power supplies are equipped with fuses for protection from overheating. The fuse responds when the temperature exceeds 40°C (104°F). Make sure power supplies have adequate cooling/ventilation. You can use screws to affix the power supply to a metal plate if necessary. Prevent build-up of excessive heat.

### **Conditions for Installation**

Before putting the system into operation for the first time, make sure there is no hazard of explosion present at the place of installation. If there is any indication that the equipment does not function properly (e.g., display remains blank, or no display backlighting) due to damage during transport, disconnect the equipment from power and notify your nearest Sartorius Service Center; see also: "Safety Inspection" in the chapter entitled "Care and Maintenance."

Take all equipment specifications into account (see "Documents for Hazardous Areas"). The explosion-protected scale must be installed and operated in accordance with acknowledged rules of engineering. Read the chapter entitled "Warnings and Safety Precautions" carefully and observe all instructions.

### **Notes on Integration into Conveyor Systems**


Any moving or rotating parts permanently attached to the EC1XS.-.....-L scale must be designed so that they cannot affect the weighing results.

For example, rotating mechanisms must be properly balanced.

The EC1XS.-.....-L must be clear on all sides during weighing so that any dirt or parts that fall will not create a connection between the weighing platform and any permanently mounted preload components. Make sure there are no cables or other objects exerting force on the load plate. Avoid generating static electricity.

### **Conditioning the Weighing System**

Moisture in the air can condense on the surface of a cold weighing instrument or other device whenever it is moved to a substantially warmer place. Condition the equipment for about 2 hours at room temperature, leaving it unplugged from AC power. Afterwards, if you keep the equipment connected to AC power, the constant positive difference in temperature between the inside of the equipment and the outside will practically rule out the effects of moisture condensation.

### **Shock Resistance**

Even though the EC1XS.-.....-L scales feature highly rugged construction, there are some limits. Avoid exposing the system to falling objects, side impact, or shocks.

### **IP Rating**

The display and control unit is rated to IP65, and the stainless steel weighing platform has IP68 protection. These ratings indicate the following levels of protection: First digit: rating 6 indicates that the equipment is dust-tight; i.e., completely resistant to penetration by solids. Second digit: rating 5 indicates resistance to penetration by water, including powerful jets of water; rating 8 indicates resistance to ingress of water during complete, continuous submersion in water to a depth of up to 10 meters (approx. 32 feet).

IP65 or IP68 protection is guaranteed only if:

- the seals on the junction box are installed correctly in accordance with industry standards, and
- the connecting cables, protective caps and cable glands were installed and connected by a qualified technician.

# Installation

## Power Supply

The scale is powered by one of the following power supplies:


- YPS05-ZD (non-hazardous area)  
or
- YPS05-XD (hazardous area)

### Note:

Observe all warnings and safety precautions. Unplug the power cord from the wall socket (mains) before performing any work on the equipment. Installation work that affects the IP65 or IP68 protection rating must be performed with extreme care. Any installation work that does not conform to the instructions in this manual will result in forfeiture of all claims under the manufacturer's warranty.

## DC Connection on the EC1XS-.....-L


- To open the display and control unit housing, remove the hexagonal nuts (7 mm) and remove the screws.


- Remove the protective cap from the bore hole on the indicator.

### Note:


Keep the seal and nut in a safe place.  
● Guide the cable with the pre-installed cable gland through the bore hole.

- Install cable gland on the housing in accordance with industry standards.
- Tighten the cable gland nut on the inside of the housing. Use a torque wrench and tighten the cable gland to 3 Nm.
- △ Check the torque at appropriate intervals to avoid over-tightening.


- Plug in the male connector on the pin terminal.
- Position the cover on the display and control unit. Check the fit of the seal between cover and housing before closing the display and control unit. Replace any defective seals. Use a torque wrench and tighten the screws to 1.6 Nm.
- Replace the hexagonal nuts (7 mm).


## Connecting the Built-in RS-232 Interface

- Open the display and control unit housing.
- Remove the protective cap.
- Connect the cable to pins 1 through 6 in accordance with industry standards (see Drawing No. 65710-740-52-A4, Sheet 2).
- Install the cable gland in accordance with industry standards.
- Position the cover on the display and control unit. Replace the hexagonal nuts (7 mm).


## Connecting the Optional RS-422 Interface

- Open the display and control unit housing.
- Remove the protective cap.
- Install the RS-422 interface card as shown in the illustration on the left.
- Connect the cable from the interface card to the main PCB (see illustration).
- Install the cable gland in accordance with industry standards.
- Connect the cable to pins 1 through 10 in accordance with industry standards (see Drawing No. 65710-740-52-A4, Sheet 4).
- Position the cover on the display and control unit. Replace the hexagonal nuts (7 mm).


## Connecting the Optional Foot Switch

- Open the display and control unit housing.
- Remove the protective cap.
- Cut the round male connector off of the model YPE05-S foot switch.
- Remove the insulation from a segment of the cable.
- Connect the cable gland (part no. YAS04CIS) in accordance with industry standards. Make sure the shield has contact to the cable gland.
- Connect the cable to the LV1 terminal strip on pins 5 and 6 in accordance with industry standards.

### Note:

Because the device is a simple switch, transposing the wires will not affect the function.

- Install the cable gland in accordance with industry standards.
- Position the cover on the display and control unit. Replace the hexagonal nuts (7 mm).

- △ Any sockets or ports on the display and control unit to which no cable is connected must be protected within hazardous areas by affixing a cap with the required IP protection rating.

# Care and Maintenance

## Service

Regular servicing by a Sartorius technician will extend the service life of your EC1XS.-.....-L scale or weighing system and help to ensure its continued weighing accuracy. Sartorius can offer you service contracts, with your choice of regular maintenance intervals ranging from 1 month to 2 years.

## Repairs

- △ Disconnect defective equipment from power immediately. Repair work must be performed by authorized Sartorius service technicians using genuine Sartorius spare parts. Any attempt by untrained persons to perform repairs may result in considerable hazards for the user.
- △ If a cable or cable gland is damaged or defective, replace the cable as a complete unit with all its connectors.

## Cleaning

- △ Always handle the equipment in keeping with its IP rating. Make sure that no liquid penetrates the equipment housing.
- Unplug the system from power before cleaning or performing any maintenance or repair work. Remove dirt and deposits from the EC1XS.-.....-L regularly.
- Wipe the equipment with a damp cloth only. Avoid generating static electricity. Devices rated to IP65 can be rinsed down with a jet of water that strikes the load plate from above.
- △ When using high-pressure cleaning equipment, aim the steam jet only at the load plate. Do not aim the jet under the scale or at the display and control unit.
- > If the water that you use to clean the equipment is too hot or too cold, the difference in temperature between the water and the equipment can cause condensation inside the equipment housing. This condensation may cause the equipment to malfunction.

## Cleaning Stainless Steel Surfaces

Clean all stainless steel parts regularly. Use a damp cloth or sponge to clean stainless steel parts on the system. You can use any commercially available household cleaning agent that is suitable for use on stainless steel. Clean stainless steel surfaces only by wiping them down. Then rinse the equipment thoroughly, making sure to remove all residues. Afterwards, wipe down stainless steel parts again using a clean, damp cloth or sponge and allow the equipment to dry. If desired, you can apply oil to the cleaned surfaces as additional protection.

- △ Do not use stainless steel cleaning agents that contain soda lye (caustic), acetic acid, hydrochloric acid, sulfuric acid or citric acid. The use of scrubbing sponges made with steel wool is not permitted. Solvents are permitted only for cleaning stainless steel parts.

## Corrosive Environment

- Remove all traces of corrosive substances on a regular basis.

## Safety Inspection

If there is any indication that safe operation of the EC1XS.-.....-L is no longer warranted; e.g.:

- if there is visible damage to the connecting cable,
- if the equipment no longer functions properly,
- if the equipment has been stored for a relatively long period under unfavorable conditions, or
- if the equipment has been subjected to rough handling during shipment,

- make sure all warnings and safety precautions are observed, and notify your nearest Sartorius Service Center or the International Technical Support Unit based in Goettingen, Germany. Maintenance and repair work may be performed only by service technicians who are authorized by Sartorius, have access to the required maintenance and repair manuals and have attended the relevant service training courses.

## Storage and Shipping Conditions

- The packaging used for shipping your Sartorius equipment is optimally designed to prevent damage during transport. It is a good idea to save the box and all parts of the packaging for future storage or shipment of the equipment. Only the original packaging provides the best protection for shipment.
- Allowable storage temperature: - 20°C to + 75°C (-4°F to +167°F)
- Allowable humidity during storage: max. 90%
- Please refer to the information under "Warnings and Safety Precautions."

## Instructions for Recycling the Packaging

If you no longer need the packaging after successful installation of the equipment, you should return it for recycling. The packaging is made from environmentally friendly materials and is a valuable source of secondary raw material. Batteries are hazardous waste and must be disposed of separately. Please deposit empty batteries in the collection boxes set up in your area for this purpose. On request, Sartorius can provide GRS boxes for collecting used batteries (GRS stands for "Gemeinsames Rücknahme System," a German organization for battery disposal\*). Contact your local waste disposal authorities if you wish to scrap the equipment. Remove batteries before scrapping the equipment. Sartorius AG in Goettingen will take back equipment and packaging for disposal in accordance with the applicable laws.\*

- \* These services are offered only within Germany. If you set up the equipment in a country other than Germany, please contact your local waste disposal authorities for information on similar services.

## Serial Number Coding

The month and year of manufacture are encoded in the serial number as follows:

YMM	x x x x x		
Y	Year		
1	2000–2006	6	2035–2041
2	2007–2013	7	2042–2048
3	2014–2020	8	2049–2055
4	2021–2027	9	2056–2062
5	2028–2034		

The first digit represents a 7-year period (see the table above); the next 2 digits represent the month. Months are numbered consecutively over the entire 7-year period, starting at 13.

2000 13-24

2001 25-36

...

Example:

113xxxxx (January 2000)

xxxxx The individual devices are numbered consecutively in the last 5 digits, starting from 00000 again at the beginning of each month.

# Declarations of Conformity

The **CE** marking may be affixed to equipment that complies with the following Directives:

## Council Directive 89/336/EEC "Electromagnetic compatibility (EMC)" Applicable European Standards:

- 1.1 Reference to 89/336/EEC: Official Journal of the European Communities, No. 2001/C105/03  
EN 61326-1 Electrical equipment for measurement, control and laboratory use  
EMC requirements  
Part 1: General requirements  
  
Defined immunity operation to interference: Industrial areas, continuous un-monitored (see Installation Instructions)  
  
Limitation of emissions: Residential areas, EC1XS.-.....-L Class B

### Note:

The operator shall be responsible for any modifications to Sartorius equipment and for any connections of cables or equipment not supplied by Sartorius. On request, Sartorius will provide information on the minimum operating specifications (in accordance with the Standards listed above).

## Council Directive 73/23/EEC "Electrical equipment designed for use within certain voltage limits"

Applicable European Standards:

- EEN 60950 Safety of information technology equipment including electrical business equipment  
EN 61010 Safety requirements for electrical equipment for measurement, control and laboratory use  
Part 1: General requirements


## 94/9/EC "Equipment and protective systems intended for use in potentially explosive atmospheres"

Applicable European Standards:


- EN 50014 General requirements  
EN 50020 Intrinsic safety

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.

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




##### Printer Code (scale operating menu)


YDP03-0CE 5-1-4 5-2-3 5-3-1 5-4-3  
YDP04IS-0CE 5-1-7 5-2-2 5-3-2 5-4-3

\*: for details on RS-232 connections, see **Sheet 3**



#### Important note!

Max. voltage  $U_m$  of any electrical apparatus connected to the YDI05-Z interface converter must not exceed 250V. Zener voltage  $U_z$  is 12.6V. Alternatively, the YDI02-ZCX barrier can be used with identical pin assignments in the 12-pin female connector; max. voltage  $U_m$  of any electrical apparatus connected to the barrier must not exceed 250V. Zener voltage  $U_z$  is 12.0V.

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
\*\*: not included with YDI05-Z

\*\*\*: see Sheet 3 for details on RS-232 connections  
A YCO01-Y may be connected to the EC1X..

Ex	Datum / Date	Name	Material	RS-422 (EX) - RS-232 (Non-EX) connection			Maßstab / Scale ---
Erstellt / Written by	12.07.05	Klausgrete			Benennung / Title <b>External connections</b>	Zeichnungs-Nr. / Drawing number <b>65710-741-50-A4</b>	
Geprüft / Reviewed by	12.07.05	Klausgrete				Blatt / Sheet <b>4</b>	
Freigabe / Released by	12.07.05	Klausgrete		Ausgabe / Revision <b>00</b>	Änderung / Alteration ---		von / of <b>4</b>

When a column is attached, the EC1XS.-.....-L can be connected to the equipotential grounding conductor either on the platform only or on the display and control unit.

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Copyright by Sartorius AG,  
Goettingen, Germany.  
All rights reserved. No part  
of this publication may  
be reprinted or translated in  
any form or by any means  
without the prior written  
permission of Sartorius AG.  
The status of the information,  
specifications and illustrations  
in this manual is indicated  
by the date given below.  
Sartorius AG reserves the  
right to make changes to  
the technology, features,  
specifications and design of  
the equipment without notice.

Status:  
September 2005, Sartorius AG,  
Goettingen, Germany