

Product data sheet

Dino-Lite AM4013MT

Dino-Lite X Pro with aluminium alloy housing (with Microtouch). The aluminium housing offers enhanced protection and durability and has the best look and feel. The AM4013MT has a built-in diffuser.



Model	AM4013MT
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	✓
LED color	White
Housing	Full aluminium alloy body
Microtouch	✓
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	✓ (line, radius, circle, 3-point circle, etc.)
Calibration	✓
Filter/Diffuser	Built-in diffuser
Compatibility	Windows XP, Vista, 7, 8 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope