

Dino-Lite Edge AM4115ZTW

The Dino-Lite Edge models provide superior image quality and a high level of flexibility due to the adaptable cap design. The built-in adjustable polarizer reduces glare and reflection on shiny objects. This model can work with 2 magnification rates while maintaining the same working distance, making this the ideal model for macro view applications with a larger field of view and a zoom function to reveal smaller details. Different caps are supplied with the microscope making the Dino-Lite Edge the perfect choice for any professional application.



Model	AM4115ZTW
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10~50x
Macro zoom	√
Interface	USB 2.0, High speed USB transmission
Nr. of LEDs	8
LED on/off switchable	√
LED color	White
Housing	Composite
Microtouch	√
Magnification lock	√
Detachable caps	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Automatic magnification reading	-
Extended Dynamic Range (EDR)	-
Extended Depth of Field (EDOF)	-
Calibration	√
Filter/Diffuser	Built-in polarizer / Diffuser cap supplied with the package
Compatibility	Windows XP, Vista, 7, 8 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Caps included in package	N3C-O- Open cap, N3C-C- Closed cap, N3C-D- Diffuser cap, N3C-L- Long cap, N3C-S- Side light cap, N3C-E Extension cap
Also available	AM4115TW without polarizer