

Reax 20 - for 4 bottles (2-32 rpm)



Overhead shaker for 4 Bottles / 2 to 32 rpm

- Set and continuously adjust the variable speed on the analog control knob from 2 to 32 rpm – ideal for waste water analysis
- Easy and fast attachments of bottles
- Other vessels with 160 - 270-mm height and max. 136-mm diameter can be used with attachments

Reax 20 - for 4 bottles (2-32 rpm) - Technical Data

Overheat protection	self-resetting
Permissible ambient conditions	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Weight	23 kg
Protection class DIN EN 60529	IP 21
Timer	-
Rotation speed range	2 - 32 rpm
Power input	280 W
Dimensions (w/d/h)	490 x 520 x 465 mm
Max. Load	30 kg
Operating Mode	-
Rotation speed setting	electronic control
Orbit / Stroke	-
Platform size	-
Accessories included	-
Angle	-
Motion	overhead